

LOINC[®] Codes for Common CDISC tests

The following table is a subset of the LOINC[®] database terms that have been identified by 3 central laboratories as the most frequently used terms in clinical trials. This is a small representation of the LOINC[®] codes and it is recommended that you use the RELMA[®] mapping tool to search the database.

Please visit the LOINC[®] web site (www.loinc.org) for more information about LOINC[®] and to download the complete set of LOINC[®] codes and the RELMA[®] mapping tool.

Explanation of fields:

Unique LOINC[®] code is a numeric code with a mod 10-check digit.

The fully specified LOINC[®] name consists of 6 fields:

1. Component (analyte) being measured
2. Property observe red (e.g. mass concentration, volume, number fraction)
3. Timing of the measurement (point in time or collected over 24 hours)
4. System (specimen)
5. Scale of measurement (quantitative, ordinal, nominal)
6. Method (where applicable)

The LOINC[®] short name is an abbreviation of the fully specified LOINC[®] name consisting of less than 30 characters. The short names are unique, but are subject to change as better algorithms for generating them are developed. It is recommended that the short name be included along with the LOINC[®] code as the 2nd part of the HL7 CE data type in HI7 messages.

LOINC [®] code	Fully Specified LOINC [®] Name						LOINC [®] short name
	Component	Property	Time Aspect	System	Scale	Method	
ALLERGY							
25638-8	EOSINOPHIL CATIONIC PROTEIN	MCNC	PT	SER	QN		ECP-mCnc
CELLMARKERS							
27011-6	CELLS.CD14	NCNC	PT	BLD	QN		CD14 Cells # Bld
8111-7	CELLS.CD14/100 CELLS	NFR	PT	BLD	QN		CD14 Cells % Bld
20402-4	CELLS.CD16+CD56+	NCNC	PT	BLD	QN		CD16+CD56+ Cells # Bld
18267-5	CELLS.CD16+CD56+/100 CELLS	NFR	PT	BLD	QN		CD16+CD56+ Cells % Bld
8116-6	CELLS.CD19	NCNC	PT	BLD	QN		CD19 Cells # Bld
8117-4	CELLS.CD19/100 CELLS	NFR	PT	BLD	QN		CD19 Cells % Bld
8122-4	CELLS.CD3	NCNC	PT	BLD	QN		CD3 Cells # Bld
8124-0	CELLS.CD3/100 CELLS	NFR	PT	BLD	QN		CD3 Cells % Bld

LOINC [®] code	Fully Specified LOINC [®] Name						LOINC [®] short name
	Component	Property	Time Aspect	System	Scale	Method	
8127-3	CELLS.CD4	NCNC	PT	BLD	QN		CD4 Cells # Bld
8128-1	CELLS.CD4/100 CELLS	NFR	PT	BLD	QN		CD4 Cells % Bld
8129-9	CELLS.CD4/CELLS.CD8	NRTO	PT	BLD	QN		CD4 Cells/CD8 Cells Bld
26982-9	CELLS.CD4+CD25+	NCNC	PT	BLD	QN		CD4+CD25+ Cells # Bld
13332-2	CELLS.CD4+CD25+/100 CELLS	NFR	PT	BLD	QN		CD4+CD25+ Cells % Bld
8137-2	CELLS.CD8	NCNC	PT	BLD	QN		CD8 Cells # Bld
8138-0	CELLS.CD8/100 CELLS	NFR	PT	BLD	QN		CD8 Cells % Bld
18349-1	CELLS.CD8+CD25+	NCNC	PT	BLD	QN		CD8+CD25+ Cells # Bld
13334-8	CELLS.CD8+CD25+/100 CELLS	NFR	PT	BLD	QN		CD8+CD25+ Cells % Bld
CHEMISTRY CHALLENGES							
14771-0	GLUCOSE^POST CFST	SCNC	PT	SER/PLAS	QN		Glucose p fast SerPl-sCnc
1558-6	GLUCOSE^POST CFST	MCNC	PT	SER/PLAS	QN		Glucose p fast SerPl-mCnc
27873-9	INSULIN^POST CFST	ACNC	PT	SER/PLAS	QN		Insulin p fast SerPl-aCnc
CHEMISTRY							
1690-7	5'-NUCLEOTIDASE	CCNC	PT	SER	QN		5NT Ser-cCnc
1715-2	ACID PHOSPHATASE	CCNC	PT	SER/PLAS	QN		ACP SerPl-cCnc
12173-1	ACID PHOSPHATASE.NON-PROSTATIC	CCNC	PT	SER	QN		ACP Non-prost Ser-cCnc
1742-6	ALANINE AMINOTRANSFERASE	CCNC	PT	SER/PLAS	QN		ALT SerPl-cCnc
1751-7	ALBUMIN	MCNC	PT	SER/PLAS	QN		Alb SerPl-mCnc
1753-3	ALBUMIN	ACNC	PT	UR	ORD		Alb Ur QI
14956-7	ALBUMIN	MRAT	24H	UR	QN	DETECTION LIMIT <= 20 MG/L	Microalbumin 24H rate Ur
14957-5	ALBUMIN	MCNC	PT	UR	QN	DETECTION LIMIT <= 20 MG/L	Microalbumin Ur Qn
30003-8	ALBUMIN	MCNC	24H	UR	QN	DETECTION LIMIT <= 20 MG/L	Microalbumin 24H Cnc Ur
2862-1	ALBUMIN	MCNC	PT	SER/PLAS	QN	ELECTROPHORESIS	Alb SerPl Elph-mCnc
14959-1	ALBUMIN/CREATININE	MCRTO	PT	UR	QN	DETECTION LIMIT <= 20 MG/L	Microalbumin/creat Ur-mRto
1759-0	ALBUMIN/GLOBULIN	MCRTO	PT	SER/PLAS	QN		A/G SerPl-mRto
13980-8	ALBUMIN/PROTEIN.TOTAL	MFR	PT	SER/PLAS	QN	ELECTROPHORESIS	Alb % SerPl Elph
13992-3	ALBUMIN/PROTEIN.TOTAL	MFR	PT	UR	QN	ELECTROPHORESIS	Alb % Ur Elph
1761-6	ALDOLASE	CCNC	PT	SER/PLAS	QN		Aldolase SerPl-cCnc
14586-2	ALDOSTERONE	SCNC	PT	SER/PLAS	QN		Aldost SerPl-sCnc
15010-2	ALDOSTERONE	SCNC	PT	UR	QN		Aldost Ur-sCnc
1763-2	ALDOSTERONE	MCNC	PT	SER/PLAS	QN		Aldost SerPl-mCnc
1764-0	ALDOSTERONE	MCNC	PT	UR	QN		Aldost Ur-mCnc
15011-0	ALDOSTERONE^SUPINE	SCNC	PT	SER/PLAS	QN		Aldost sup SerPl-sCnc
1767-3	ALDOSTERONE^SUPINE	MCNC	PT	SER/PLAS	QN		Aldost sup SerPl-mCnc
15012-8	ALDOSTERONE^UPRIGHT	SCNC	PT	SER/PLAS	QN		Aldost upr SerPl-sCnc
1768-1	ALDOSTERONE^UPRIGHT	MCNC	PT	SER/PLAS	QN		Aldost upr SerPl-mCnc
6768-6	ALKALINE PHOSPHATASE	CCNC	PT	SER/PLAS	QN		ALP SerPl-cCnc
1777-2	ALKALINE PHOSPHATASE.BONE	CCNC	PT	SER/PLAS	QN		ALP Bone SerPl-cCnc
15013-6	ALKALINE PHOSPHATASE.BONE/ALKALINE PHOSPHATASE.TOTAL	CFR	PT	SER/PLAS	QN		ALP Bone % SerPl
1778-0	ALKALINE PHOSPHATASE.INTESTINAL	CCNC	PT	SER/PLAS	QN		ALP Intest SerPl-cCnc
15014-4	ALKALINE PHOSPHATASE.INTESTINAL/ALKALINE PHOSPHATASE.TOTAL	CFR	PT	SER/PLAS	QN		ALP Intest % SerPl
1779-8	ALKALINE PHOSPHATASE.LIVER	CCNC	PT	SER/PLAS	QN		ALP Liver SerPl-cCnc

LOINC [®] code	Fully Specified LOINC [®] Name						LOINC [®] short name
	Component	Property	Time Aspect	System	Scale	Method	
15015-1	ALKALINE PHOSPHATASE.LIVER/ALKALINE PHOSPHATASE.TOTAL	CFR	PT	SER/PLAS	QN		ALP Liver % SerPl
1825-9	ALPHA 1 ANTITRYPSIN	MCNC	PT	SER/PLAS	QN		A1AT SerPl-mCnc
2865-4	ALPHA 1 GLOBULIN	MCNC	PT	SER/PLAS	QN	ELECTROPHORESIS	A1 Globulin SerPl Elph-mCnc
13978-2	ALPHA 1 GLOBULIN/PROTEIN.TOTAL	MFR	PT	SER/PLAS	QN	ELECTROPHORESIS	A1 Globulin % SerPl Elph
13990-7	ALPHA 1 GLOBULIN/PROTEIN.TOTAL	MFR	PT	UR	QN	ELECTROPHORESIS	A1 Globulin % Ur Elph
2868-8	ALPHA 2 GLOBULIN	MCNC	PT	SER/PLAS	QN	ELECTROPHORESIS	A2 Globulin SerPl Elph-mCnc
13981-6	ALPHA 2 GLOBULIN/PROTEIN.TOTAL	MFR	PT	SER/PLAS	QN	ELECTROPHORESIS	A2 Globulin % SerPl Elph
13993-1	ALPHA 2 GLOBULIN/PROTEIN.TOTAL	MFR	PT	UR	QN	ELECTROPHORESIS	A2 Globulin % Ur Elph
14590-4	ALPHA TOCOPHEROL	SCNC	PT	SER	QN		A-Tocopherol Vit E Ser-sCnc
2685-6	ALPHA-1-ACID GLYCOPROTEIN	MCNC	PT	SER/PLAS	QN		A1Glycoprot SerPl-mCnc
1834-1	ALPHA-1-FETOPROTEIN	MCNC	PT	SER	QN		AFP Ser-mCnc
2644-3	ALPHA-N-ACETYLGLUCOSAMINIDASE	CCNC	PT	UR	QN		A-NAG Ur-cCnc
16362-6	AMMONIA	SCNC	PT	PLAS	QN		Ammonia Plas-sCnc
22763-7	AMMONIA	MCNC	PT	PLAS	QN		Ammonia Plas-mCnc
1798-8	AMYLASE	CCNC	PT	SER/PLAS	QN		Amylase SerPl-cCnc
1805-1	AMYLASE.PANCREATIC	CCNC	PT	SER/PLAS	QN		Amylase P SerPl-cCnc
15057-3	ANDROSTANOLONE	SCNC	PT	SER/PLAS	QN		Androstanolone SerPl-sCnc
1848-1	ANDROSTANOLONE	MCNC	PT	SER/PLAS	QN		Androstanolone SerPl-mCnc
14603-5	ANDROSTENEDIONE	SCNC	PT	SER/PLAS	QN		Androst SerPl-sCnc
1854-9	ANDROSTENEDIONE	MCNC	PT	SER/PLAS	QN		Androst SerPl-mCnc
1869-7	APOLIPOPROTEIN A-I	MCNC	PT	SER/PLAS	QN		Apo A-I SerPl-mCnc
1870-5	APOLIPOPROTEIN A-II	MCNC	PT	SER/PLAS	QN		Apo A-II SerPl-mCnc
1884-6	APOLIPOPROTEIN B	MCNC	PT	SER/PLAS	QN		Apo B SerPl-mCnc
1877-0	APOLIPOPROTEIN C-III	MCNC	PT	SER/PLAS	QN		Apo C-III SerPl-mCnc
1886-1	APOLIPOPROTEIN E	MCNC	PT	SER/PLAS	QN		Apo E SerPl-mCnc
1920-8	ASPARTATE AMINOTRANSFERASE	CCNC	PT	SER/PLAS	QN		AST SerPl-cCnc
2871-2	BETA GLOBULIN	MCNC	PT	SER/PLAS	QN	ELECTROPHORESIS	B-Globulin SerPl Elph-mCnc
13982-4	BETA GLOBULIN/PROTEIN.TOTAL	MFR	PT	SER/PLAS	QN	ELECTROPHORESIS	B-Globulin % SerPl Elph
13994-9	BETA GLOBULIN/PROTEIN.TOTAL	MFR	PT	UR	QN	ELECTROPHORESIS	B-Globulin % Ur Elph
29512-1	BETA HYDROXYBUTYRATE	MCNC	PT	SER/PLAS	QN		B-OH-Butyr SerPl-mCnc
6873-4	BETA HYDROXYBUTYRATE	SCNC	PT	SER/PLAS	QN		B-OH-Butyr SerPl-sCnc
1951-3	BETA-2-MICROGLOBULIN	MCNC	PT	CSF	QN		B2 Microglob CSF-mCnc
1952-1	BETA-2-MICROGLOBULIN	MCNC	PT	SER	QN		B2 Microglob Ser-mCnc
1953-9	BETA-2-MICROGLOBULIN	MCNC	PT	UR	QN		B2 Microglob Ur-mCnc
1959-6	BICARBONATE	SCNC	PT	BLD	QN		HCO3 Bld-sCnc
1963-8	BICARBONATE	SCNC	PT	SER	QN		HCO3 Ser-sCnc
14628-2	BILE ACID	SCNC	PT	SER	QN		Bile Ac Ser-sCnc
14631-6	BILIRUBIN	SCNC	PT	SER/PLAS	QN		Bilirub SerPl-sCnc
1974-5	BILIRUBIN	MCNC	PT	FLU	QN		Bilirub Fld-mCnc
1975-2	BILIRUBIN	MCNC	PT	SER/PLAS	QN		Bilirub SerPl-mCnc
1977-8	BILIRUBIN	ACNC	PT	UR	ORD		Bilirub Ur QI
29767-1	BILIRUBIN	SCNC	PT	FLU	QN		Bilirub Fld-sCnc
14629-0	BILIRUBIN.GLUCURONIDATED+BILIRUBIN.ALBUMIN BOUND	SCNC	PT	SER/PLAS	QN		Bilirub Direct SerPl-sCnc
1968-7	BILIRUBIN.GLUCURONIDATED+BILIRUBIN.ALBUMIN BOUND	MCNC	PT	SER/PLAS	QN		Bilirub Direct SerPl-mCnc
14630-8	BILIRUBIN.NON-GLUCURONIDATED	SCNC	PT	SER/PLAS	QN		Bilirub Indirect SerPl-sCnc
1971-1	BILIRUBIN.NON-GLUCURONIDATED	MCNC	PT	SER/PLAS	QN		Bilirub Indirect SerPl-mCnc
14633-2	C PEPTIDE	SCNC	PT	SER/PLAS	QN		C Peptide SerPl-sCnc

LOINC [®] code	Fully Specified LOINC [®] Name						LOINC [®] short name
	Component	Property	Time Aspect	System	Scale	Method	
1986-9	C PEPTIDE	MCNC	PT	SER/PLAS	QN		C Peptide SerPl-mCnc
1988-5	C REACTIVE PROTEIN	MCNC	PT	SER/PLAS	QN		CRP SerPl-mCnc
30522-7	C REACTIVE PROTEIN	MCNC	PT	SER/PLAS	QN	HIGH SENSITIVITY	CRP SerPl High Sens-mCnc
14635-7	CALCIDIOL	SCNC	PT	SER/PLAS	QN		Vit D25 SerPl-sCnc
1989-3	CALCIDIOL	MCNC	PT	SER/PLAS	QN		Vit D25 SerPl-mCnc
14637-3	CALCIUM	SRAT	24H	UR	QN		Calcium 24H Ur-sRate
17861-6	CALCIUM	MCNC	PT	SER/PLAS	QN		Calcium SerPl-mCnc
17862-4	CALCIUM	MCNC	PT	UR	QN		Calcium Ur-mCnc
18488-7	CALCIUM	MCNC	24H	UR	QN		Calcium 24H Ur-mCnc
2000-8	CALCIUM	SCNC	PT	SER/PLAS	QN		Calcium SerPl-sCnc
2004-0	CALCIUM	SCNC	PT	UR	QN		Calcium Ur-sCnc
25362-5	CALCIUM	SCNC	24H	UR	QN		Calcium 24H Ur-sCnc
1995-0	CALCIUM.IONIZED	SCNC	PT	SER/PLAS	QN		Ca-I SerPl-sCnc
18281-6	CALCIUM^^CORRECTED FOR TOTAL PROTEIN	SCNC	PT	SER	QN		Calcium TP cor Ser-sCnc
24108-3	CANCER AG 19-9	ACNC	PT	SER/PLAS	QN		Cancer Ag19-9 SerPl-aCnc
2028-9	CARBON DIOXIDE	SCNC	PT	SER/PLAS	QN		CO2 SerPl-sCnc
20563-3	CARBOXYHEMOGLOBIN/HEMOGLOBIN.TOTAL	MFR	PT	BLD	QN		COHgb % Bld
2039-6	CARCINOEMBRYONIC AG	MCNC	PT	SER/PLAS	QN		CEA SerPl-mCnc
14287-7	CARNITINE	SCNC	PT	UR	QN		Carnitine Ur-sCnc
14288-5	CARNITINE	SCNC	PT	SER/PLAS	QN		Carnitine SerPl-sCnc
24448-3	CARNITINE	SRAT	24H	UR	QN		Carnitine 24H Ur-sRate
14285-1	CARNITINE.FREE	SCNC	PT	UR	QN		Carnitine Free Ur-sCnc
14286-9	CARNITINE.FREE	SCNC	PT	SER/PLAS	QN		Carnitine Free SerPl-sCnc
2064-4	CERULOPLASMIN	MCNC	PT	SER	QN		Ceruloplasmin Ser-mCnc
2075-0	CHLORIDE	SCNC	PT	SER/PLAS	QN		Chloride SerPl-sCnc
2078-4	CHLORIDE	SCNC	PT	UR	QN		Chloride Ur-sCnc
21194-6	CHLORIDE	SCNC	24H	UR	QN		Chloride 24H Ur-sCnc
14647-2	CHOLESTEROL	SCNC	PT	SER/PLAS	QN		Cholest SerPl-sCnc
2093-3	CHOLESTEROL	MCNC	PT	SER/PLAS	QN		Cholest SerPl-mCnc
14646-4	CHOLESTEROL.IN HDL	SCNC	PT	SER/PLAS	QN		HDLc SerPl-sCnc
2085-9	CHOLESTEROL.IN HDL	MCNC	PT	SER/PLAS	QN		HDLc SerPl-mCnc
26017-4	CHOLESTEROL.IN HDL 3	SCNC	PT	SER/PLAS	QN		HDL3c SerPl-sCnc
9833-5	CHOLESTEROL.IN HDL 3	MCNC	PT	SER/PLAS	QN		HDL3c SerPl-mCnc
2089-1	CHOLESTEROL.IN LDL	MCNC	PT	SER/PLAS	QN		LDLc SerPl-mCnc
22748-8	CHOLESTEROL.IN LDL	SCNC	PT	SER/PLAS	QN		LDLc SerPl-sCnc
18262-6	CHOLESTEROL.IN LDL	MCNC	PT	SER/PLAS	QN	DIRECT ASSAY	LDLc SerPl Direct Assay-mCnc
2091-7	CHOLESTEROL.IN VLDL	MCNC	PT	SER/PLAS	QN		VLDLc SerPl-mCnc
25371-6	CHOLESTEROL.IN VLDL	SCNC	PT	SER/PLAS	QN		VLDLc SerPl-sCnc
2110-5	CHORIOGONADOTROPIN.BETA SUBUNIT	ACNC	PT	SER	ORD		B-HCG Ser QI
21198-7	CHORIOGONADOTROPIN.BETA SUBUNIT	ACNC	PT	SER	QN		B-HCG Ser-aCnc
6687-8	CITRATE	MRAT	24H	UR	QN		Citrate 24H Ur-mRate
14685-2	COBALAMINS	SCNC	PT	SER	QN		Vit B12 Ser-sCnc
2132-9	COBALAMINS	MCNC	PT	SER	QN		Vit B12 Ser-mCnc
21215-9	COLLAGEN CROSSLINKED N-TELOPEPTIDE	SCNC	PT	SER	QN		Collagen NTx Ser-sCnc
27939-8	COLLAGEN CROSSLINKED N-TELOPEPTIDE	SCNC	PT	UR	QN		Collagen NTx Ur-sCnc
14115-0	COLLAGEN CROSSLINKED N-TELOPEPTIDE/CREATININE	SCRTO	PT	UR	QN		Collagen NTx/creat Ur-sRto
14674-6	CORTICOTROPIN	SCNC	PT	PLAS	QN		ACTH Plas-sCnc
2141-0	CORTICOTROPIN	MCNC	PT	PLAS	QN		ACTH Plas-mCnc
14675-3	CORTISOL	SCNC	PT	SER/PLAS	QN		Cortis SerPl-sCnc
2143-6	CORTISOL	MCNC	PT	SER/PLAS	QN		Cortis SerPl-mCnc

LOINC® code	Fully Specified LOINC® Name						LOINC® short name
	Component	Property	Time Aspect	System	Scale	Method	
25883-0	CORTISOL.FREE	SCNC	24H	UR	QN	HPLC	Cortis Free 24H Ur HPLC-sCnc
14679-5	CORTISOL^MORNING	SCNC	PT	SER/PLAS	QN		Cortis AM SerPl-sCnc
2157-6	CREATINE KINASE	CCNC	PT	SER/PLAS	QN		CK SerPl-cCnc
2152-7	CREATINE KINASE.BB	CCNC	PT	SER/PLAS	QN	ELECTROPHORESIS	CK BB SerPl Elph-cCnc
9642-0	CREATINE KINASE.BB/CREATINE KINASE.TOTAL	CFR	PT	SER/PLAS	QN		CK BB % SerPl
13969-1	CREATINE KINASE.MB	MCNC	PT	SER/PLAS	QN		CK MB SerPl-mCnc
6773-6	CREATINE KINASE.MB	CCNC	PT	SER/PLAS	QN	EIA	CK MB SerPl EIA-cCnc
2154-3	CREATINE KINASE.MB	CCNC	PT	SER/PLAS	QN	ELECTROPHORESIS	CK MB SerPl Elph-cCnc
20569-0	CREATINE KINASE.MB/CREATINE KINASE.TOTAL	CFR	PT	SER/PLAS	QN		CK MB % SerPl
2155-0	CREATINE KINASE.MM	CCNC	PT	SER/PLAS	QN	ELECTROPHORESIS	CK MM SerPl Elph-cCnc
9643-8	CREATINE KINASE.MM/CREATINE KINASE.TOTAL	CFR	PT	SER/PLAS	QN		CK MM % SerPl
14682-9	CREATININE	SCNC	PT	SER/PLAS	QN		Creat SerPl-sCnc
14683-7	CREATININE	SCNC	PT	UR	QN		Creat Ur-sCnc
14684-5	CREATININE	SRAT	24H	UR	QN		Creat 24H Ur-sRate
15051-6	CREATININE	SCNC	PT	DIAF	QN		Creat Diaf-sCnc
20624-3	CREATININE	MCNC	24H	UR	QN		Creat 24H Ur-mCnc
2160-0	CREATININE	MCNC	PT	SER/PLAS	QN		Creat SerPl-mCnc
2161-8	CREATININE	MCNC	PT	UR	QN		Creat Ur-mCnc
2162-6	CREATININE	MRAT	24H	UR	QN		Creat 24H Ur-mRate
25886-3	CREATININE	SCNC	24H	UR	QN		Creat 24H Ur-sCnc
5919-6	CREATININE	MCNC	PT	DIAF	QN		Creat Diaf-mCnc
15054-0	DEHYDROEPIANDROSTERONE	SCNC	PT	SER/PLAS	QN		DHEA SerPl-sCnc
2193-1	DEHYDROEPIANDROSTERONE	MCNC	PT	SER/PLAS	QN		DHEA SerPl-mCnc
14688-6	DEHYDROEPIANDROSTERONE SULFATE	SCNC	PT	SER/PLAS	QN		DHEA-S SerPl-sCnc
2191-5	DEHYDROEPIANDROSTERONE SULFATE	MCNC	PT	SER/PLAS	QN		DHEA-S SerPl-mCnc
27424-1	DEOXYPYRIDINOLINE	SCNC	PT	UR	QN		DPD Ur-sCnc
25095-1	DEOXYPYRIDINOLINE/CREATININE	SCRTO	PT	UR	QN		DPD/creat Ur-sRto
15061-5	ERYTHROPOIETIN	ACNC	PT	SER/PLAS	QN		EPO SerPl-aCnc
14715-7	ESTRADIOL	SCNC	PT	SER/PLAS	QN		Estradiol SerPl-sCnc
2243-4	ESTRADIOL	MCNC	PT	SER/PLAS	QN		Estradiol SerPl-mCnc
2254-1	ESTROGEN	MCNC	PT	SER/PLAS	QN		Estrogen SerPl-mCnc
22663-9	ESTRONE	SCNC	PT	SER/PLAS	QN		Estrone SerPl-sCnc
15066-4	FATTY ACIDS.NONESTERIFIED	SCNC	PT	SER/PLAS	QN		NEFA SerPl-sCnc
14723-1	FERRITIN	SCNC	PT	SER	QN		Ferritin Ser-sCnc
2276-4	FERRITIN	MCNC	PT	SER	QN		Ferritin Ser-mCnc
14732-2	FOLATE	SCNC	PT	SER	QN		Folate Ser-sCnc
2284-8	FOLATE	MCNC	PT	SER	QN		Folate Ser-mCnc
15067-2	FOLLITROPIN	ACNC	PT	SER/PLAS	QN		FSH SerPl-aCnc
15069-8	FRUCTOSAMINE	SCNC	PT	SER/PLAS	QN		Fructosamine SerPl-sCnc
2874-6	GAMMA GLOBULIN	MCNC	PT	SER/PLAS	QN	ELECTROPHORESIS	G-Globulin SerPl Elph-mCnc
13983-2	GAMMA GLOBULIN/PROTEIN.TOTAL	MFR	PT	SER/PLAS	QN	ELECTROPHORESIS	G-Globulin % SerPl Elph
13995-6	GAMMA GLOBULIN/PROTEIN.TOTAL	MFR	PT	UR	QN	ELECTROPHORESIS	G-Globulin % Ur Elph
2324-2	GAMMA GLUTAMYL TRANSFERASE	CCNC	PT	SER/PLAS	QN		GGT SerPl-cCnc
2333-3	GASTRIN	MCNC	PT	SER/PLAS	QN		Gastrin SerPl-mCnc
2336-6	GLOBULIN	MCNC	PT	SER	QN		Globulin Ser-mCnc
14749-6	GLUCOSE	SCNC	PT	SER/PLAS	QN		Glucose SerPl-sCnc
15076-3	GLUCOSE	SCNC	PT	UR	QN		Glucose Ur-sCnc
2345-7	GLUCOSE	MCNC	PT	SER/PLAS	QN		Glucose SerPl-mCnc
2349-9	GLUCOSE	ACNC	PT	UR	ORD		Glucose Ur QI
2350-7	GLUCOSE	MCNC	PT	UR	QN		Glucose Ur-mCnc
2415-8	HISTAMINE	MCNC	PT	BLD	QN		Histamine Bld-mCnc

LOINC [®] code	Fully Specified LOINC [®] Name						LOINC [®] short name
	Component	Property	Time Aspect	System	Scale	Method	
13965-9	HOMOCYSTEINE	SCNC	PT	SER/PLAS	QN		Homocysteine SerPl-sCnc
2458-8	IGA	MCNC	PT	SER	QN		IgA Ser-mCnc
19113-0	IGE	ACNC	PT	SER	QN		IgE Ser-aCnc
2465-3	IGG	MCNC	PT	SER	QN		IgG Ser-mCnc
2472-9	IGM	MCNC	PT	SER	QN		IgM Ser-mCnc
20448-7	INSULIN	ACNC	PT	SER/PLAS	QN		Insulin SerPl-aCnc
25447-4	INSULIN.FREE	SCNC	PT	SER/PLAS	QN		Insulin Free SerPl-sCnc
2483-6	INSULIN-LIKE GROWTH FACTOR BINDING PROTEIN 3	MCNC	PT	SER/PLAS	QN		IGF BP3 SerPl-mCnc
2484-4	INSULIN-LIKE GROWTH FACTOR-I	MCNC	PT	SER/PLAS	QN		IGF-I SerPl-mCnc
29775-4	INSULIN-LIKE GROWTH FACTOR-I	SCNC	PT	SER/PLAS	QN		IGF-I SerPl-sCnc
13440-3	INTERPRETATION	IMP	PT	UR	NOM	IMMUNOFIXATION	Interpretation Ur IFE-Imp
14798-3	IRON	SCNC	PT	SER/PLAS	QN		Iron SerPl-sCnc
2498-4	IRON	MCNC	PT	SER/PLAS	QN		Iron SerPl-mCnc
25937-4	IRON	SCNC	24H	UR	QN		Iron 24H Ur-sCnc
30030-1	IRON	MCNC	24H	UR	QN		Iron 24H Ur-mCnc
14800-7	IRON BINDING CAPACITY	SCNC	PT	SER/PLAS	QN		TIBC SerPl-sCnc
2500-7	IRON BINDING CAPACITY	MCNC	PT	SER/PLAS	QN		TIBC SerPl-mCnc
22753-8	IRON BINDING CAPACITY.UNSATURATED	SCNC	PT	SER/PLAS	QN		UIBC SerPl-sCnc
2501-5	IRON BINDING CAPACITY.UNSATURATED	MCNC	PT	SER/PLAS	QN		UIBC SerPl-mCnc
2514-8	KETONES	ACNC	PT	UR	ORD	TEST STRIP	Ketones Ur Qi Strip
14118-4	LACTATE	MCNC	PT	SER/PLAS	QN		Lactate SerPl-mCnc
2524-7	LACTATE	SCNC	PT	SER/PLAS	QN		Lactate SerPl-sCnc
2532-0	LACTATE DEHYDROGENASE	CCNC	PT	SER/PLAS	QN		LDH SerPl-cCnc
2534-6	LACTATE DEHYDROGENASE	CCNC	PT	UR	QN		LDH Ur-cCnc
2537-9	LACTATE DEHYDROGENASE 1	CCNC	PT	SER/PLAS	QN		LDH1 SerPl-cCnc
2540-3	LACTATE DEHYDROGENASE 2	CCNC	PT	SER/PLAS	QN		LDH2 SerPl-cCnc
2543-7	LACTATE DEHYDROGENASE 3	CCNC	PT	SER/PLAS	QN		LDH3 SerPl-cCnc
2546-0	LACTATE DEHYDROGENASE 4	CCNC	PT	SER/PLAS	QN		LDH4 SerPl-cCnc
2549-4	LACTATE DEHYDROGENASE 5	CCNC	PT	SER/PLAS	QN		LDH5 SerPl-cCnc
21365-2	LEPTIN	MCNC	PT	SER/PLAS	QN		Leptin SerPl-mCnc
10835-7	LIPOPROTEIN (LITTLE A)	MCNC	PT	SER/PLAS	QN		LPa SerPl-mCnc
10501-5	LUTROPIN	ACNC	PT	SER/PLAS	QN		LH SerPl-aCnc
19123-9	MAGNESIUM	MCNC	PT	SER/PLAS	QN		Magnesium SerPl-mCnc
19124-7	MAGNESIUM	MCNC	PT	UR	QN		Magnesium Ur-mCnc
25954-9	MAGNESIUM	SCNC	24H	UR	QN		Magnesium 24H Ur-sCnc
2598-1	MAGNESIUM	SCNC	PT	UR	QN		Magnesium Ur-sCnc
2601-3	MAGNESIUM	SCNC	PT	SER/PLAS	QN		Magnesium SerPl-sCnc
2639-3	MYOGLOBIN	MCNC	PT	SER/PLAS	QN		Myoglobin SerPl-mCnc
2641-9	MYOGLOBIN	MCNC	PT	UR	QN		Myoglobin Ur-mCnc
30934-4	NATRIURETIC PEPTIDE.B	MCNC	PT	SER/PLAS	QN		BNP SerPl-mCnc
9740-2	NEOPTERIN	SCNC	PT	SER/PLAS	QN		Neopterin SerPl-sCnc
2695-5	OSMOLALITY	OSMOL	PT	UR	QN		Osmolality Ur Qn
2692-2	OSMOLALITY	OSMOL	PT	SER/PLAS	QN		Osmolality SerPl Qn
2694-8	OSMOLALITY	OSMOL	24H	UR	QN		Osmolality 24H Ur Qn
15084-7	OSTEOCALCIN	SCNC	PT	SER	QN		Osteocalcin Ser-sCnc
2697-1	OSTEOCALCIN	MCNC	PT	SER	QN		Osteocalcin Ser-mCnc
2713-6	OXYGEN SATURATION.CALC FROM OXYGEN PARTIAL PRESSURE	SFR	PT	BLD	QN		O2 Satn from pO2 % Bld
14866-8	PARATHYRIN.INTACT	SCNC	PT	SER/PLAS	QN		PTH-Intact SerPl-sCnc
2731-8	PARATHYRIN.INTACT	MCNC	PT	SER/PLAS	QN		PTH-Intact SerPl-mCnc
27378-9	PH	SCNC	24H	UR	QN		pH 24H Ur-sCnc

LOINC [®] code	Fully Specified LOINC [®] Name						LOINC [®] short name
	Component	Property	Time Aspect	System	Scale	Method	
2753-2	PH	SCNC	PT	SER/PLAS	QN		pH SerPI-sCnc
2756-5	PH	SCNC	PT	UR	QN		pH Ur-sCnc
13539-2	PHOSPHATE	SCNC	PT	UR	QN		Phosphate Ur-sCnc
14879-1	PHOSPHATE	SCNC	PT	SER/PLAS	QN		Phosphate SerPI-sCnc
14881-7	PHOSPHATE	SRAT	24H	UR	QN		Phosphate 24H Ur-sRate
21458-5	PHOSPHATE	MCNC	24H	UR	QN		Phosphate 24H Ur-mCnc
25973-9	PHOSPHATE	SCNC	24H	UR	QN		Phosphate 24H Ur-sCnc
2775-5	PHOSPHATE	MCNC	PT	DIAF	QN		Phosphate Diaf-mCnc
2777-1	PHOSPHATE	MCNC	PT	SER/PLAS	QN		Phosphate SerPI-mCnc
2778-9	PHOSPHATE	MCNC	PT	UR	QN		Phosphate Ur-mCnc
2779-7	PHOSPHATE	MRAT	24H	UR	QN		Phosphate 24H Ur-mRate
21476-7	POTASSIUM	SCNC	24H	UR	QN		Potassium 24H Ur-sCnc
2823-3	POTASSIUM	SCNC	PT	SER/PLAS	QN		Potassium SerPI-sCnc
2828-2	POTASSIUM	SCNC	PT	UR	QN		Potassium Ur-sCnc
2829-0	POTASSIUM	SRAT	24H	UR	QN		Potassium 24H Ur-sRate
14890-8	PROGESTERONE	SCNC	PT	SER/PLAS	QN		Progest SerPI-sCnc
2839-9	PROGESTERONE	MCNC	PT	SER/PLAS	QN		Progest SerPI-mCnc
2842-3	PROLACTIN	MCNC	PT	SER/PLAS	QN		Prolactin SerPI-mCnc
2857-1	PROSTATE SPECIFIC AG	MCNC	PT	SER/PLAS	QN		PSA SerPI-mCnc
21482-5	PROTEIN	MCNC	24H	UR	QN		Prot 24H Ur-mCnc
2885-2	PROTEIN	MCNC	PT	SER/PLAS	QN		Prot SerPI-mCnc
2889-4	PROTEIN	MRAT	24H	UR	QN		Prot 24H Ur-mRate
2890-2	PROTEIN/CREATININE	MCRTO	PT	UR	QN		Prot/creat Ur-mRto
25977-0	PYRIDINOLINE	SCNC	24H	UR	QN		PYD 24H Ur-sCnc
25129-8	PYRIDINOLINE/CREATININE	SCRTO	PT	UR	QN		PYD/creat Ur-sRto
13967-5	SEX HORMONE BINDING GLOBULIN	SCNC	PT	SER	QN		SHBG Ser-sCnc
21525-1	SODIUM	SCNC	24H	UR	QN		Sodium 24H Ur-sCnc
2951-2	SODIUM	SCNC	PT	SER/PLAS	QN		Sodium SerPI-sCnc
2955-3	SODIUM	SCNC	PT	UR	QN		Sodium Ur-sCnc
2956-1	SODIUM	SRAT	24H	UR	QN		Sodium 24H Ur-sRate
16257-8	SODIUM URATE CRYSTALS	ACNC	PT	CALC	ORD	INFRARED SPECTROSCOPY	Na Urate Cry Calc QI IR
2963-7	SOMATOTROPIN	MCNC	PT	SER/PLAS	QN		GH SerPI-mCnc
2965-2	SPECIFIC GRAVITY	RDEN	PT	UR	QN		Sp Gr Ur Qn
14913-8	TESTOSTERONE	SCNC	PT	SER/PLAS	QN		Testost SerPI-sCnc
2986-8	TESTOSTERONE	MCNC	PT	SER/PLAS	QN		Testost SerPI-mCnc
14914-6	TESTOSTERONE.FREE	SCNC	PT	SER/PLAS	QN		Testost Free SerPI-sCnc
2991-8	TESTOSTERONE.FREE	MCNC	PT	SER/PLAS	QN		Testost Free SerPI-mCnc
11061-9	THROMBOXANE BETA 2	MRAT	24H	UR	QN		TXB2 24H Ur-mRate
3016-3	THYROTROPIN	ACNC	PT	SER/PLAS	QN		TSH SerPI-aCnc
14921-1	THYROXINE	SCNC	PT	SER/PLAS	QN		T4 SerPI-sCnc
27980-2	THYROXINE BINDING GLOBULIN	SCNC	PT	SER/PLAS	QN		T4BG SerPI-sCnc
3021-3	THYROXINE BINDING GLOBULIN	MCNC	PT	SER/PLAS	QN		T4BG SerPI-mCnc
32215-6	THYROXINE FREE INDEX	ACNC	PT	SER/PLAS	QN		FTI SerPI-aCnc
14920-3	THYROXINE.FREE	SCNC	PT	SER/PLAS	QN		Free T4 SerPI-sCnc
3024-7	THYROXINE.FREE	MCNC	PT	SER/PLAS	QN		Free T4 SerPI-mCnc
22674-6	TRANSFERRIN	SCNC	PT	SER/PLAS	QN		Transferrin SerPI-sCnc
3034-6	TRANSFERRIN	MCNC	PT	SER/PLAS	QN		Transferrin SerPI-mCnc
30248-9	TRANSFERRIN RECEPTOR.SOLUBLE	MCNC	PT	SER/PLAS	QN		sTfR SerPI-mCnc
3040-3	TRIACYLGLYCEROL LIPASE	CCNC	PT	SER/PLAS	QN		Lipase SerPI-cCnc
14927-8	TRIGLYCERIDE	SCNC	PT	SER/PLAS	QN		Trigl SerPI-sCnc
2571-8	TRIGLYCERIDE	MCNC	PT	SER/PLAS	QN		Trigl SerPI-mCnc

LOINC® code	Fully Specified LOINC® Name						LOINC® short name
	Component	Property	Time Aspect	System	Scale	Method	
14930-2	TRIIODOTHYRONINE	SCNC	PT	SER/PLAS	QN		T3 SerPl-sCnc
3053-6	TRIIODOTHYRONINE	MCNC	PT	SER/PLAS	QN		T3 SerPl-mCnc
3050-2	TRIIODOTHYRONINE RESIN UPTAKE	NFR	PT	SER/PLAS	QN		T3RU % SerPl
14928-6	TRIIODOTHYRONINE.FREE	SCNC	PT	SER/PLAS	QN		T3Free SerPl-sCnc
3051-0	TRIIODOTHYRONINE.FREE	MCNC	PT	SER/PLAS	QN		T3Free SerPl-mCnc
14929-4	TRIIODOTHYRONINE.REVERSE	SCNC	PT	SER/PLAS	QN		T3Reverse SerPl-sCnc
3052-8	TRIIODOTHYRONINE.REVERSE	MCNC	PT	SER/PLAS	QN		T3Reverse SerPl-mCnc
10839-9	TROPONIN I.CARDIAC	MCNC	PT	SER/PLAS	QN		Troponin I SerPl-mCnc
6598-7	TROPONIN T.CARDIAC	MCNC	PT	SER/PLAS	QN		Troponin T SerPl-mCnc
21582-2	TRYPTASE	MCNC	PT	SER/PLAS	QN		Tryptase SerPl-mCnc
13900-6	TUMOR NECROSIS FACTOR.ALPHA	MCNC	PT	SER	QN		A-Tumor Necrosis Fact Ser-mCnc
14933-6	URATE	SCNC	PT	SER/PLAS	QN		Urate SerPl-sCnc
14934-4	URATE	SCNC	PT	UR	QN		Urate Ur-sCnc
21587-1	URATE	MCNC	24H	UR	QN		Urate 24H Ur-mCnc
25997-8	URATE	SCNC	24H	UR	QN		Urate 24H Ur-sCnc
3084-1	URATE	MCNC	PT	SER/PLAS	QN		Urate SerPl-mCnc
3086-6	URATE	MCNC	PT	UR	QN		Urate Ur-mCnc
22664-7	UREA	SCNC	PT	SER/PLAS	QN		Urea SerPl-sCnc
14936-9	UREA NITROGEN	SCNC	PT	DIAFP	QN		Urea Nit DiafP-sCnc
14937-7	UREA NITROGEN	SCNC	PT	SER/PLAS	QN		BUN SerPl-sCnc
14938-5	UREA NITROGEN	SCNC	PT	UR	QN		UUN Ur-sCnc
3094-0	UREA NITROGEN	MCNC	PT	SER/PLAS	QN		BUN SerPl-mCnc
3095-7	UREA NITROGEN	MCNC	PT	UR	QN		UUN Ur-mCnc
5918-8	UREA NITROGEN	MCNC	PT	DIAF	QN		Urea Nit Diaf-mCnc
3097-3	UREA NITROGEN/CREATININE	MCRTO	PT	SER/PLAS	QN		BUN/creat SerPl-mRto
13658-0	UROBILINOGEN	ACNC	PT	UR	ORD		Urobilinogen Ur QI
3107-0	UROBILINOGEN	MCNC	PT	UR	QN		Urobilinogen Ur-mCnc
COAGULATION							
27811-9	ANTITHROMBIN ACTUAL/NORMAL	RLCCNC	PT	PPP	QN	ENZY	AT III Act/Nor PPP Chro
27812-7	ANTITHROMBIN AG ACTUAL/NORMAL	RLMCNC	PT	PPP	QN	IMM	AT III Ag Act/Nor PPP Imm
8065-5	CARDIOLIPIN AB.IGG	ACNC	PT	SER	QN		Cardiolipin IgG Ser-aCnc
8067-1	CARDIOLIPIN AB.IGM	ACNC	PT	SER	QN		Cardiolipin IgM Ser-aCnc
3218-5	COAGULATION FACTOR X ACTIVITY ACTUAL/NORMAL	RLTM	PT	PPP	QN	COAG	FX Act/Nor PPP Qn
5904-8	COAGULATION RUSSELL VIPER VENOM INDUCED	TIME	PT	PPP	QN	COAG	Coag RVV Ind PPP Qn
13488-2	COAGULATION SURFACE INDUCED	TIME	PT	PPP^CONTROL	QN	COAG	aPTT PPP Cont Qn
14979-9	COAGULATION SURFACE INDUCED	TIME	PT	PPP	QN	COAG	aPTT PPP Qn
3173-2	COAGULATION SURFACE INDUCED	TIME	PT	BLD	QN	COAG	aPTT Bld Qn
5946-9	COAGULATION SURFACE INDUCED.FACTOR SUBSTITUTION^IMMEDIATELY AFTER ADDITION OF NORMAL PLASMA	TIME	PT	PPP	QN	COAG	aPTT imm NP PPP Qn
3243-3	COAGULATION THROMBIN INDUCED	TIME	PT	PPP	QN	COAG	TT PPP Qn
5901-4	COAGULATION TISSUE FACTOR INDUCED	TIME	PT	PPP^CONTROL	QN	COAG	PT PPP Cont Qn
5902-2	COAGULATION TISSUE FACTOR INDUCED	TIME	PT	PPP	QN	COAG	PT PPP Qn
6301-6	COAGULATION TISSUE FACTOR INDUCED.INR	RLTM	PT	PPP	QN	COAG	INR PPP Qn
30240-6	FIBRIN D-DIMER	MCNC	PT	PPP	QN		D Dimer PPP-mCnc
3255-7	FIBRINOGEN	MCNC	PT	PPP	QN	COAG	Fibrinogen PPP-mCnc
22758-7	PLASMINOGEN ACTIVATOR INHIBITOR 1 AG	MCNC	PT	PPP	QN	IMM	PAI1 Ag PPP Imm-mCnc
24377-4	PLATELET AGGREGATION.ADENOSINE DIPHOSPHATE INDUCED	ACNC	PT	PRP	ORD		IPA ADP PRP QI

LOINC [®] code	Fully Specified LOINC [®] Name						LOINC [®] short name
	Component	Property	Time Aspect	System	Scale	Method	
24379-0	PLATELET AGGREGATION.COLLAGEN INDUCED	ACNC	PT	PRP	ORD		IPA Coll PRP QI
14182-0	THROMBIN ANTITHROMBIN COMPLEX AG	MCNC	PT	PPP	QN	IMM	TAT Ag PPP Imm-mCnc
27816-8	VON WILLEBRAND FACTOR AG ACTUAL/NORMAL	RLMCNC	PT	PPP	QN	IMM	vWf Ag Act/Nor PPP Imm
CYTOLOGY							
19767-3	CYTOLOGIST	ID	PT	CVX/VAG	NOM	CYTO STAIN	Cytologist Cvx/Vag Cyto
19774-9	CYTOLOGY STUDY COMMENT	IMP	PT	CVX/VAG	NAR	CYTO STAIN	Cytology Cmnt Cvx/Vag Cyto-IMP
19762-4	GENERAL CATEGORIES	IMP	PT	CVX/VAG	NOM	CYTO STAIN	Gen Categ Cvx/Vag Cyto- Imp
19769-9	PATHOLOGIST	ID	PT	CVX/VAG	NOM	CYTO STAIN	Pathologist Cvx/Vag Cyto
19773-1	RECOMMENDED FOLLOWUP	PRID	PT	CVX/VAG	NOM	CYTO STAIN	Recom F/U Cvx/Vag Cyto
19763-2	SPECIMEN SOURCE	PRID	PT	CVX/VAG	NOM	CYTO STAIN	Spec Source Cvx/Vag Cyto
19764-0	STATEMENT OF ADEQUACY	IMP	PT	CVX/VAG	NOM	CYTO STAIN	Stat of Adq Cvx/Vag Cyto- Imp
DRUG LEVELS & TOXICOLOGY							
8150-5	AMPHETAMINES	MCNC	PT	UR	QN		Amphets Ur-mCnc
16369-1	AMPHETAMINES	ACNC	PT	UR	ORD	CONFIRM	Amphets Ur QI Cfm
19261-7	AMPHETAMINES	ACNC	PT	UR	ORD	SCREEN	Amphets Ur QI Scn
9426-8	BARBITURATES	MCNC	PT	UR	QN		Barbiturates Ur-mCnc
16429-3	BARBITURATES	ACNC	PT	UR	ORD	CONFIRM	Barbiturates Ur QI Cfm
19270-8	BARBITURATES	ACNC	PT	UR	ORD	SCREEN	Barbiturates Ur QI Scn
3390-2	BENZODIAZEPINES	ACNC	PT	UR	ORD		Benzodiaz Ur QI
9428-4	BENZODIAZEPINES	MCNC	PT	UR	QN		Benzodiaz Ur-mCnc
16195-0	BENZODIAZEPINES	ACNC	PT	UR	ORD	CONFIRM	Benzodiaz Ur QI Cfm
14316-4	BENZODIAZEPINES	ACNC	PT	UR	ORD	SCREEN	Benzodiaz Ur QI Scn
14315-6	BENZOYLECGONINE	ACNC	PT	UR	ORD	CONFIRM	BZE Ur QI Cfm
14314-9	BENZOYLECGONINE	ACNC	PT	UR	ORD	SCREEN	BZE Ur QI Scn
19289-8	CANNABINOIDS	ACNC	PT	UR	ORD	CONFIRM	Cannabinoids Ur QI Cfm
18282-4	CANNABINOIDS	ACNC	PT	UR	ORD	SCREEN	Cannabinoids Ur QI Scn
14639-9	CARBAMAZEPINE	SCNC	PT	SER/PLAS	QN		CBZ SerPl-sCnc
3432-2	CARBAMAZEPINE	MCNC	PT	SER/PLAS	QN		CBZ SerPl-mCnc
3398-5	COCAINE	MCNC	PT	UR	QN		Cocaine Ur-mCnc
10365-5	COTININE	MCNC	PT	SER/PLAS	QN		Cotinine SerPl-mCnc
12293-7	COTININE	ACNC	PT	UR	ORD		Cotinine Ur QI
9372-4	COTININE	MCNC	PT	XXX	QN		Cotinine XXX-mCnc
3530-3	DELTA-9-TETRAHYDROCANNABINOL	MCNC	PT	UR	QN		D9THC Ur-mCnc
19419-1	DEXTROAMPHETAMINE	ACNC	PT	UR	ORD	SCREEN	D-amphet Ur QI Scn
10535-3	DIGOXIN	MCNC	PT	SER/PLAS	QN		Digoxin SerPl-mCnc
14698-5	DIGOXIN	SCNC	PT	SER/PLAS	QN		Digoxin SerPl-sCnc
14719-9	ETHANOL	SCNC	PT	SER/PLAS	QN		Ethanol SerPl-sCnc
5639-0	ETHANOL	ACNC	PT	BLD	ORD		Ethanol Bld QI
5640-8	ETHANOL	MCNC	PT	BLD	ORD		Ethanol Bld-mCnc
5643-2	ETHANOL	MCNC	PT	SER/PLAS	QN		Ethanol SerPl-mCnc
5644-0	ETHANOL	ACNC	PT	UR	ORD		Ethanol Ur QI
14334-7	LITHIUM	SCNC	PT	SER/PLAS	QN		Lithium SerPl-sCnc
3719-2	LITHIUM	MCNC	PT	SER/PLAS	QN		Lithium SerPl-mCnc
19528-9	LYSERGATE DIETHYLAMIDE	ACNC	PT	UR	ORD	SCREEN	LSD Ur QI Scn
3774-7	METHADONE	MCNC	PT	UR	QN		Methadone Ur-mCnc
16199-2	METHADONE	ACNC	PT	UR	ORD	CONFIRM	Methadone Ur QI Cfm
19550-3	METHADONE	ACNC	PT	UR	ORD	SCREEN	Methadone Ur QI Scn

LOINC® code	Fully Specified LOINC® Name						LOINC® short name
	Component	Property	Time Aspect	System	Scale	Method	
19554-5	METHAMPHETAMINE	ACNC	PT	UR	ORD	SCREEN	Methamphet Ur Ql Scn
3787-9	METHAQUALONE	MCNC	PT	UR	QN		Methaqualone Ur-mCnc
18389-7	METHAQUALONE	ACNC	PT	UR	ORD	CONFIRM	Methaqualone Ur Ql Cfm
19558-6	METHAQUALONE	ACNC	PT	UR	ORD	SCREEN	Methaqualone Ur Ql Scn
8220-6	OPIATES	MCNC	PT	UR	QN		Opiates Ur-mCnc
18390-5	OPIATES	ACNC	PT	UR	ORD	CONFIRM	Opiates Ur Ql Cfm
19295-5	OPIATES	ACNC	PT	UR	ORD	SCREEN	Opiates Ur Ql Scn
3937-0	PHENCYCLIDINE	MCNC	PT	UR	QN		PCP Ur-mCnc
18392-1	PHENCYCLIDINE	ACNC	PT	UR	ORD	CONFIRM	PCP Ur Ql Cfm
19659-2	PHENCYCLIDINE	ACNC	PT	UR	ORD	SCREEN	PCP Ur Ql Scn
3954-5	PHENOTHIAZINE	MCNC	PT	UR	QN		Phenothiaz Ur-mCnc
16221-4	PHENOTHIAZINE	ACNC	PT	UR	ORD	CONFIRM	Phenothiaz Ur Ql Cfm
19670-9	PHENOTHIAZINE	ACNC	PT	UR	ORD	SCREEN	Phenothiaz Ur Ql Scn
3545-1	PROPOXYPHENE	MCNC	PT	UR	QN		Propoxyph Ur-mCnc
16200-8	PROPOXYPHENE	ACNC	PT	UR	ORD	CONFIRM	Propoxyph Ur Ql Cfm
19429-0	PROPOXYPHENE	ACNC	PT	UR	ORD	SCREEN	Propoxyph Ur Ql Scn
14915-3	THEOPHYLLINE	SCNC	PT	SER/PLAS	QN		Theophylline SerPl-sCnc
4049-3	THEOPHYLLINE	MCNC	PT	SER/PLAS	QN		Theophylline SerPl-mCnc
22752-0	TOBRAMYCIN^TROUGH	SCNC	PT	SER/PLAS	QN		Tobramycin Tr SerPl-sCnc
4059-2	TOBRAMYCIN^TROUGH	MCNC	PT	SER/PLAS	QN		Tobramycin Tr SerPl-mCnc
16181-0	TRICYCLIC ANTIDEPRESSANTS	MCNC	PT	UR	QN	CONFIRM	Tricyclics Ur Cfm-mCnc
19315-1	TRICYCLIC ANTIDEPRESSANTS	ACNC	PT	UR	ORD	CONFIRM	Tricyclics Ur Ql Cfm
19312-8	TRICYCLIC ANTIDEPRESSANTS	ACNC	PT	UR	ORD	SCREEN	Tricyclics Ur Ql Scn
HISTORY							
8665-2	DATE LAST MENSTRUAL PERIOD	TMSTP	PT	^PATIENT	QN	REPORTED	
HEMATOLOGY & CELL COUNTS							
7789-1	ACANTHOCYTES	ACNC	PT	BLD	ORD	MICROSCOPY.LIGHT	Acanthocytes Bld Ql Smear
702-1	ANISOCYTOSIS	ACNC	PT	BLD	ORD	MICROSCOPY.LIGHT	Anisocytosis Bld Ql Smear
11281-3	AUER RODS	ACNC	PT	BLD	ORD	MICROSCOPY.LIGHT	Auer Bodies Bld Ql Smear
703-9	BASOPHILIC STIPPLING	ACNC	PT	BLD	ORD	MICROSCOPY.LIGHT	Baso Stipl Bld Ql Smear
26444-0	BASOPHILS	NCNC	PT	BLD	QN		Basophils # Bld
30180-4	BASOPHILS/100 LEUKOCYTES	NFR	PT	BLD	QN		Basophils % Bld
10371-3	BITE CELLS	ACNC	PT	BLD	ORD	MICROSCOPY.LIGHT	Bite Cells Bld Ql Smear
30376-8	BLASTS	NCNC	PT	BLD	QN		Blasts # Bld
26446-5	BLASTS/100 LEUKOCYTES	NFR	PT	BLD	QN		Blasts % Bld
7790-9	BURR CELLS	ACNC	PT	BLD	ORD	MICROSCOPY.LIGHT	Burr Cells Bld Ql Smear
11280-5	CABOT RINGS	ACNC	PT	BLD	ORD	MICROSCOPY.LIGHT	Cabot Rings Bld Ql Smear
4485-9	COMPLEMENT C3	MCNC	PT	SER/PLAS	QN		C3 SerPl-mCnc
4498-2	COMPLEMENT C4	MCNC	PT	SER/PLAS	QN		C4 SerPl-mCnc
4532-8	COMPLEMENT TOTAL HEMOLYTIC CH50	ACNC	PT	SER/PLAS	QN		CH50 SerPl-aCnc
7791-7	DACRYOCYTES	ACNC	PT	BLD	ORD	MICROSCOPY.LIGHT	Dacryocytes Bld Ql Smear
7792-5	DOHLE BODY	ACNC	PT	BLD	ORD	MICROSCOPY.LIGHT	Dohle Bod Bld Ql Smear
11274-8	ELLIPTOCYTES	ACNC	PT	BLD	ORD	MICROSCOPY.LIGHT	Elliptocytes Bld Ql Smear
26449-9	EOSINOPHILS	NCNC	PT	BLD	QN		Eosinophil # Bld
26450-7	EOSINOPHILS/100 LEUKOCYTES	NFR	PT	BLD	QN		Eosinophil % Bld
30385-9	ERYTHROCYTE DISTRIBUTION WIDTH	RATIO	PT	RBC	QN		RDW RBC-Rto
28539-5	ERYTHROCYTE MEAN CORPUSCULAR HEMOGLOBIN	ENTMASS	PT	RBC	QN		MCH RBC Qn
28540-3	ERYTHROCYTE MEAN CORPUSCULAR HEMOGLOBIN CONCENTRATION	MCNC	PT	RBC	QN		MCHC RBC-mCnc
4537-7	ERYTHROCYTE SEDIMENTATION RATE	VEL	PT	BLD	QN	WESTERGREN	ESR Bld Qn Westrgm

LOINC [®] code	Fully Specified LOINC [®] Name						LOINC [®] short name
	Component	Property	Time Aspect	System	Scale	Method	
26453-1	ERYTHROCYTES	NCNC	PT	BLD	QN		RBC # Bld
26461-4	ERYTHROCYTES.NUCLEATED/100 ERYTHROCYTES	NFR	PT	BLD	QN		nRBC/RBC % Bld
19048-8	ERYTHROCYTES.NUCLEATED/100 LEUKOCYTES	RATIO	PT	BLD	QN		nRBC/WBC Bld-Rto
30397-4	HAIRY CELLS	NCNC	PT	BLD	QN		Hairy Cells # Bld
24106-7	HAIRY CELLS/100 LEUKOCYTES	NFR	PT	BLD	QN	MANUAL COUNT	Hairy Cells % Bld Manual
4542-7	HAPTOGLOBIN	MCNC	PT	SER	QN		Haptoglob Ser-mCnc
716-1	HEINZ BODIES	ACNC	PT	BLD	ORD	MICROSCOPY.LIGHT	Heinz Bod Bld Ql Smear
10374-7	HELMET CELLS	ACNC	PT	BLD	ORD	MICROSCOPY.LIGHT	Helmet Cells Bld Ql Smear
20570-8	HEMATOCRIT	VFR	PT	BLD	QN		Hct % Bld
4545-0	HEMATOCRIT	VFR	PT	BLD	QN	SPUN	Hct % Bld Spun
14134-1	HEMOGLOBIN	MCNC	PT	SER	QN		Hgb Ser-mCnc
25433-4	HEMOGLOBIN	SCNC	PT	PLAS	QN		Hgb Plas-sCnc
718-7	HEMOGLOBIN	MCNC	PT	BLD	QN		Hgb Bld-mCnc
20572-4	HEMOGLOBIN A/HEMOGLOBIN.TOTAL	MFR	PT	BLD	QN	ELECTROPHORESIS	Hgb A % Bld Elph
4547-6	HEMOGLOBIN A1/HEMOGLOBIN.TOTAL	MFR	PT	BLD	QN		Hgb A1 % Bld
4548-4	HEMOGLOBIN A1C/HEMOGLOBIN.TOTAL	MFR	PT	BLD	QN		Hgb A1c % Bld
17856-6	HEMOGLOBIN A1C/HEMOGLOBIN.TOTAL	MFR	PT	BLD	QN	HPLC	Hgb A1c % Bld HPLC
4562-5	HEMOGLOBIN C/HEMOGLOBIN.TOTAL	MFR	PT	BLD	QN	ELECTROPHORESIS PH 8.9	Hgb C % Bld Elph Alk
4625-0	HEMOGLOBIN S/HEMOGLOBIN.TOTAL	MFR	PT	BLD	QN		Hgb S % Bld
7793-3	HOWELL-JOLLY BODIES	ACNC	PT	BLD	ORD	MICROSCOPY.LIGHT	Howell-Jolly Bod Bld Ql Smear
30400-6	HYPOCHROMIA	ACNC	PT	BLD	ORD		Hypochromia Bld Ql
27415-9	INTERFERON.GAMMA	MCNC	PT	SER/PLAS	QN		G-IFN SerPl-mCnc
13629-1	INTERLEUKIN 1 BETA	MCNC	PT	SER/PLAS	QN		I1Beta SerPl-mCnc
26848-2	INTERLEUKIN 10	MCNC	PT	SER/PLAS	QN		I10 SerPl-mCnc
27161-9	INTERLEUKIN 4	MCNC	PT	SER/PLAS	QN		I4 SerPl-mCnc
26881-3	INTERLEUKIN 6	MCNC	PT	SER/PLAS	QN		I6 SerPl-mCnc
26462-2	LARGE UNSTAINED CELLS	NCNC	PT	BLD	QN		LUC # Bld
26463-0	LARGE UNSTAINED CELLS/100 LEUKOCYTES	NFR	PT	BLD	QN		LUC % Bld
26464-8	LEUKOCYTES	NCNC	PT	BLD	QN		WBC # Bld
30405-5	LEUKOCYTES	NCNC	PT	UR	QN		WBC # Ur
26474-7	LYMPHOCYTES	NCNC	PT	BLD	QN		Lymphocytes # Bld
26477-0	LYMPHOCYTES.ATYPICAL	NCNC	PT	BLD	QN		Atypical Lymphs # Bld
13046-8	LYMPHOCYTES.ATYPICAL/100 LEUKOCYTES	NFR	PT	BLD	QN		Atypical Lymphs % Bld
26478-8	LYMPHOCYTES/100 LEUKOCYTES	NFR	PT	BLD	QN		Lymphocytes % Bld
30422-0	LYMPHOMA CELLS	NCNC	PT	BLD	QN		Lymphoma Cells # Bld
30423-8	LYMPHOMA CELLS/100 LEUKOCYTES	NFR	PT	BLD	QN		Lymphoma Cells % Bld
30424-6	MACROCYTES	ACNC	PT	BLD	ORD		Macrocytes Bld Ql
30428-7	MEAN CORPUSCULAR VOLUME	ENTVOL	PT	RBC	QN		MCV RBC Qn
30433-7	METAMYELOCYTES	NCNC	PT	BLD	QN		Metamyelocytes # Bld
28541-1	METAMYELOCYTES/100 LEUKOCYTES	NFR	PT	BLD	QN		Metamyelocytes % Bld
30434-5	MICROCYTES	ACNC	PT	BLD	ORD		Microcytes Bld Ql
26484-6	MONOCYTES	NCNC	PT	BLD	QN		Monocytes # Bld
26485-3	MONOCYTES/100 LEUKOCYTES	NFR	PT	BLD	QN		Monocytes % Bld
18314-5	MORPHOLOGY	IMP	PT	BLD	NAR		Morphology Bld-Imp
30444-4	MYELOBLASTS	NCNC	PT	BLD	QN		Myeloblasts # Bld
30445-1	MYELOBLASTS/100 LEUKOCYTES	NFR	PT	BLD	QN		Myeloblasts % Bld
30446-9	MYELOCYTES	NCNC	PT	BLD	QN		Myelocytes # Bld
26498-6	MYELOCYTES/100 LEUKOCYTES	NFR	PT	BLD	QN		Myelocytes % Bld
26499-4	NEUTROPHILS	NCNC	PT	BLD	QN		Neutrophils # Bld
26507-4	NEUTROPHILS.BAND FORM	NCNC	PT	BLD	QN		Neuts Band # Bld

LOINC® code	Fully Specified LOINC® Name						LOINC® short name
	Component	Property	Time Aspect	System	Scale	Method	
26508-2	NEUTROPHILS.BAND FORM/100 LEUKOCYTES	NFR	PT	BLD	QN		Neuts Band % Bld
765-8	NEUTROPHILS.HYPERSEGMENTED	ACNC	PT	BLD	ORD	MICROSCOPY.LIGHT	Neuts Hyperseg Bld Ql Smear
30451-9	NEUTROPHILS.SEGMENTED	NCNC	PT	BLD	QN		Neuts Seg # Bld
26505-8	NEUTROPHILS.SEGMENTED/100 LEUKOCYTES	NFR	PT	BLD	QN		Neuts Seg % Bld
18319-4	NEUTROPHILS.VACUOLATED	ACNC	PT	BLD	ORD	MICROSCOPY.LIGHT	Neuts Vac Bld Ql Smear
26511-6	NEUTROPHILS/100 LEUKOCYTES	NFR	PT	BLD	QN		Neutrophils % Bld
774-0	OVALOCYTES	ACNC	PT	BLD	ORD	MICROSCOPY.LIGHT	Ovalocytes Bld Ql Smear
7795-8	PAPPENHEIMER BODIES	ACNC	PT	BLD	ORD	MICROSCOPY.LIGHT	Pappenheimer Bod Bld Ql Smear
18311-1	PELGER HUET CELLS	ACNC	PT	BLD	ORD	MICROSCOPY.LIGHT	Pelger Huet Cells Bld Ql Smear
30458-4	PLASMA CELLS	NCNC	PT	BLD	QN		Plasma Cells # Bld
11117-9	PLASMA CELLS.IMMATURE/100 CELLS	NFR	PT	MAR	QN	MICROSCOPY	Imm Plasma Cells % Mar Micro
13047-6	PLASMA CELLS/100 LEUKOCYTES	NFR	PT	BLD	QN		Plasma Cells % Bld
28542-9	PLATELET MEAN VOLUME	ENTVOL	PT	BLD	QN		PMV Bld Qn
18312-9	PLATELET SATELLITISM	ACNC	PT	BLD	ORD	MICROSCOPY.LIGHT	Platelet Satel Bld Ql Smear
26515-7	PLATELETS	NCNC	PT	BLD	QN		Platelet # Bld
5908-9	PLATELETS.GIANT	ACNC	PT	BLD	ORD	MICROSCOPY.LIGHT	Giant Platelets Bld Ql Smear
32146-3	PLATELETS.LARGE	ACNC	PT	BLD	ORD	MICROSCOPY.LIGHT	Large Platelets Bld Ql Smear
779-9	POIKILOCYTOSIS	ACNC	PT	BLD	ORD	MICROSCOPY.LIGHT	Poikilocytosis Bld Ql Smear
10378-8	POLYCHROMASIA	ACNC	PT	BLD	ORD	MICROSCOPY.LIGHT	Polychromasia Bld Ql Smear
30464-2	PROLYMPHOCYTES	NCNC	PT	BLD	QN		Prolymphocytes # Bld
30465-9	PROLYMPHOCYTES/100 LEUKOCYTES	NFR	PT	BLD	QN		Prolymphocytes % Bld
30466-7	PROMONOCYTES/100 LEUKOCYTES	NFR	PT	BLD	QN		Promonocytes % Bld
26523-1	PROMYELOCYTES	NCNC	PT	BLD	QN		Promyelocytes # Bld
26524-9	PROMYELOCYTES/100 LEUKOCYTES	NFR	PT	BLD	QN		Promyelocytes % Bld
14196-0	RETICULOCYTES	NCNC	PT	RBC	QN		Retics # RBC
17849-1	RETICULOCYTES/100 ERYTHROCYTES	NFR	PT	RBC	QN	AUTOMATED COUNT	Retics/100 RBC % Auto
7797-4	ROULEAUX	ACNC	PT	BLD	ORD	MICROSCOPY.LIGHT	Rouleaux Bld Ql Smear
800-3	SCHISTOCYTES	ACNC	PT	BLD	ORD	MICROSCOPY.LIGHT	Schistocytes Bld Ql Smear
801-1	SICKLE CELLS	ACNC	PT	BLD	ORD	MICROSCOPY.LIGHT	Sickle Cells Bld Ql Smear
7798-2	SMUDGE CELLS	ACNC	PT	BLD	ORD	MICROSCOPY.LIGHT	Smudge Cells Bld Ql Smear
802-9	SPHEROCYTES	ACNC	PT	BLD	ORD	MICROSCOPY.LIGHT	Spherocytes Bld Ql Smear
10380-4	STOMATOCYTES	ACNC	PT	BLD	ORD	MICROSCOPY.LIGHT	Stomatocytes Bld Ql Smear
10381-2	TARGET CELLS	ACNC	PT	BLD	ORD	MICROSCOPY.LIGHT	Targets Bld Ql Smear
803-7	TOXIC GRANULES	ACNC	PT	BLD	ORD	MICROSCOPY.LIGHT	Toxic Granules Bld Ql Smear
MICROBIOLOGY							
22203-4	CLOSTRIDIUM TETANI AB.IGG	ACNC	PT	SER	QN		C tetani IgG Ser-aCnc
13227-4	CORYNEBACTERIUM DIPHTHERIAE AB.IGG	ACNC	PT	SER	QN		C diphtheriae IgG Ser-aCnc
13949-3	CYTOMEGALOVIRUS AB.IGG	ACNC	PT	SER	ORD	EIA	CMV IgG Ser Ql EIA
20575-7	HEPATITIS A VIRUS AB	ACNC	PT	SER	ORD		Hep A Ab Ser Ql
22314-9	HEPATITIS A VIRUS AB.IGM	ACNC	PT	SER	ORD		Hep A IgM Ser Ql
16933-4	HEPATITIS B VIRUS CORE AB	ACNC	PT	SER	ORD		Hep Bc Ab Ser Ql
22319-8	HEPATITIS B VIRUS CORE AB.IGM	ACNC	PT	SER	QN		Hep Bc IgM Ser-aCnc
29900-8	HEPATITIS B VIRUS DNA	MCNC	PT	SER	QN	PROBE.AMP.SIG	Hep B DNA Ser bDNA-mCnc
29615-2	HEPATITIS B VIRUS DNA	NCNC	PT	SER/PLAS	QN	PROBE.AMP.TAR	Hep B Viral Load SerPl PCR
22320-6	HEPATITIS B VIRUS LITTLE E AB	ACNC	PT	SER	ORD		Hep Be Ab Ser Ql
13954-3	HEPATITIS B VIRUS LITTLE E AG	ACNC	PT	SER	ORD	EIA	Hep Be Ag Ser Ql EIA

LOINC® code	Fully Specified LOINC® Name						LOINC® short name
	Component	Property	Time Aspect	System	Scale	Method	
22322-2	HEPATITIS B VIRUS SURFACE AB	ACNC	PT	SER	ORD		Hep Bs Ab Ser Ql
5195-3	HEPATITIS B VIRUS SURFACE AG	ACNC	PT	SER	ORD		Hep Bs Ag Ser Ql
5196-1	HEPATITIS B VIRUS SURFACE AG	ACNC	PT	SER	ORD	EIA	Hep Bs Ag Ser Ql EIA
7905-3	HEPATITIS B VIRUS SURFACE AG	ACNC	PT	SER	ORD	NEUT	Hep Bs Ag Ser Ql Nt
16128-1	HEPATITIS C VIRUS AB	ACNC	PT	SER	ORD		Hep C Ab Ser Ql
11259-9	HEPATITIS C VIRUS RNA	ACNC	PT	SER/PLAS	ORD	PROBE.AMP.TAR	Hep C RNA SerPl Ql PCR
20416-4	HEPATITIS C VIRUS RNA	NCNC	PT	SER/PLAS	QN	PROBE.AMP.TAR	Hep C Viral Load SerPl PCR
13248-0	HEPATITIS D VIRUS AB	ACNC	PT	SER	ORD		Hep D Ab Ser Ql
7909-5	HERPES SIMPLEX VIRUS 1 AB.IGG	ACNC	PT	SER	QN		HSV1 IgG Ser-aCnc
7912-9	HERPES SIMPLEX VIRUS 2 AB.IGG	ACNC	PT	SER	QN		HSV2 IgG Ser-aCnc
19106-4	HERPES SIMPLEX VIRUS AB.IGG	ACNC	PT	SER	ORD		HSV IgG Ser Ql
7917-8	HIV 1 AB	ACNC	PT	SER	ORD		HIV1 Ab Ser Ql
20447-9	HIV 1 RNA	NCNC	PT	SER/PLAS	QN	PROBE.AMP.TAR	HIV1 Viral Load SerPl PCR
7918-6	HIV 1+2 AB	ACNC	PT	SER	ORD		HIV1+2 Ab Ser Ql
7919-4	HIV 2 AB	ACNC	PT	SER	ORD		HIV2 Ab Ser Ql
6447-7	LEGIONELLA PNEUMOPHILA AG	ACNC	PT	UR	ORD	EIA	L pneumo Ag Ur Ql EIA
21406-4	MYCOPLASMA PNEUMONIAE AB.IGM	ACNC	PT	SER	ORD		M pneumo IgM Ser Ql
29675-6	PARVOVIRUS B19 AB.IGG	ACNC	PT	SER	ORD		B19V IgG Ser Ql
10710-2	PLASMODIUM SP IDENTIFIED	PRID	PT	BLD	NOM	THIN FILM	Plasmodium Bld Thin Film
20507-0	REAGIN AB	ACNC	PT	SER	ORD	RAPID TEST	RPR Ser Ql
5370-2	STREPTOLYSIN O AB	ACNC	PT	SER	QN		ASO Ab Ser-aCnc
6561-5	TREPONEMA PALLIDUM AB.IGG	ACNC	PT	SER	ORD		T pallidum IgG Ser Ql
6562-3	TREPONEMA PALLIDUM AB.IGM	ACNC	PT	SER	ORD		T pallidum IgM Ser Ql
SEROLOGY							
11013-0	DNA DOUBLE STRAND AB	TITR	PT	SER	QN		dsDNA Ab Titr Ser
5130-0	DNA DOUBLE STRAND AB	ACNC	PT	SER	QN		dsDNA Ab Ser-aCnc
12277-0	DNA DOUBLE STRAND AB	ACNC	PT	SER	ORD	FARR	dsDNA Ab Ser Ql Farr
5131-8	DNA DOUBLE STRAND AB	ACNC	PT	SER	ORD	IF	dsDNA Ab Ser Ql IF
8061-4	NUCLEAR AB	ACNC	PT	SER	ORD		ANA Ser Ql
14611-8	NUCLEAR AB PATTERN	IMP	PT	SER	NOM		ANA Pat Ser-Imp
11572-5	RHEUMATOID FACTOR	ACNC	PT	SER	QN		RF Ser-aCnc
SPECIMEN							
5778-6	COLOR	TYPE	PT	UR	NOM		Color Ur
4691-2	VISCOSITY	VISC	PT	PLAS	QN		Viscosity Plas Qn
URINALYSIS							
5766-1	AMMONIUM URATE CRYSTALS	ACNC	PT	URNS	ORD	MICROSCOPY.LIGHT	Amm Urate Cry UrnS Ql
25145-4	BACTERIA	ACNC	PT	URNS	ORD	MICROSCOPY.LIGHT	Bacteria UrnS Ql
5771-1	BILIRUBIN CRYSTALS	ACNC	PT	URNS	ORD	MICROSCOPY.LIGHT	Bilirub Cry UrnS Ql
5773-7	CALCIUM CARBONATE CRYSTALS	ACNC	PT	URNS	ORD	MICROSCOPY.LIGHT	Ca Carbonate Cry UrnS Ql
5774-5	CALCIUM OXALATE CRYSTALS	ACNC	PT	URNS	ORD	MICROSCOPY.LIGHT	CaOx Cry UrnS Ql
5775-2	CALCIUM PHOSPHATE CRYSTALS	ACNC	PT	URNS	ORD	MICROSCOPY.LIGHT	Ca Phos Cry UrnS Ql
24124-0	CASTS	ACNC	PT	URNS	ORD	MICROSCOPY.LIGHT	Casts UrnS Ql
9842-6	CASTS	NARIC	PT	URNS	QN	MICROSCOPY.LIGHT.LPF	Casts UrnS Qn LPF
5777-8	CHOLESTEROL CRYSTALS	ACNC	PT	URNS	ORD	MICROSCOPY.LIGHT	Cholest Cry UrnS Ql
32167-9	CLARITY	TYPE	PT	UR	NOM		Clarity Ur
5782-8	CRYSTALS	PRID	PT	URNS	NOM	MICROSCOPY.LIGHT	Crystals UrnS
5784-4	CYSTINE CRYSTALS	ACNC	PT	URNS	ORD	MICROSCOPY.LIGHT	Cystine Cry UrnS Ql
25157-9	EPITHELIAL CASTS	ACNC	PT	URNS	ORD	MICROSCOPY.LIGHT	Epith Casts UrnS Ql

LOINC [®] code	Fully Specified LOINC [®] Name						LOINC [®] short name
	Component	Property	Time Aspect	System	Scale	Method	
5786-9	EPITHELIAL CASTS	NARIC	PT	URNS	QN	MICROSCOPY.LIGHT.LPF	Epith Casts UrnS Qn LPF
20453-7	EPITHELIAL CELLS	ACNC	PT	URNS	ORD	MICROSCOPY.LIGHT	Epi Cells UrnS Ql
5787-7	EPITHELIAL CELLS	NARIC	PT	URNS	QN	MICROSCOPY.LIGHT.HPF	Epi Cells UrnS Qn HPF
12248-1	EPITHELIAL CELLS.RENAL	ACNC	PT	URNS	ORD	MICROSCOPY.LIGHT	Renal Epi Cells UrnS Ql
5807-3	ERYTHROCYTE CASTS	NARIC	PT	URNS	QN	MICROSCOPY.LIGHT.LPF	RBC Casts UrnS Qn LPF
14290-1	ERYTHROCYTES	ACNC	PT	URNS	QN	MANUAL COUNT	RBC UrnS Manual-aCnc
13945-1	ERYTHROCYTES	NARIC	PT	URNS	QN	MICROSCOPY.LIGHT.HPF	RBC UrnS Qn HPF
5808-1	ERYTHROCYTES	NCNC	PT	URNS	QN	MICROSCOPY.LIGHT.HPF	RBC # UrnS HPF
20409-9	ERYTHROCYTES	NCNC	PT	UR	QN	TEST STRIP	RBC # Ur Strip
5788-5	FAT BODIES.OVAL	NARIC	PT	URNS	QN	MICROSCOPY.LIGHT.HPF	Oval Fat Bod UrnS Qn HPF
25159-5	FATTY CASTS	ACNC	PT	URNS	ORD	MICROSCOPY.LIGHT	Fatty Casts UrnS Ql
5789-3	FATTY CASTS	NARIC	PT	URNS	QN	MICROSCOPY.LIGHT.LPF	Fatty Casts UrnS Qn LPF
25160-3	GRANULAR CASTS	ACNC	PT	URNS	ORD	MICROSCOPY.LIGHT	Gran Casts UrnS Ql
5793-5	GRANULAR CASTS	NARIC	PT	URNS	QN	MICROSCOPY.LIGHT.LPF	Gran Casts UrnS Qn LPF
5794-3	HEMOGLOBIN	ACNC	PT	UR	ORD	TEST STRIP	Hgb Ur Ql Strip
5795-0	HIPPURATE CRYSTALS	ACNC	PT	URNS	ORD	MICROSCOPY.LIGHT	Hippurate Cry UrnS Ql
25162-9	HYALINE CASTS	ACNC	PT	URNS	ORD	MICROSCOPY.LIGHT	Hyaline Casts UrnS Ql
5796-8	HYALINE CASTS	NARIC	PT	URNS	QN	MICROSCOPY.LIGHT.LPF	Hyaline Casts UrnS Qn LPF
5798-4	LEUCINE CRYSTALS	ACNC	PT	URNS	ORD	MICROSCOPY.LIGHT	Leucine Cry UrnS Ql
5799-2	LEUKOCYTE ESTERASE	ACNC	PT	UR	ORD	TEST STRIP	WBC Est Ur Ql Strip
5821-4	LEUKOCYTES	NARIC	PT	URNS	QN	MICROSCOPY.LIGHT.HPF	WBC UrnS Qn HPF
12235-8	MICROSCOPIC OBSERVATION	PRID	PT	URNS	NOM	MICROSCOPY.LIGHT	Micro UrnS
8247-9	MUCUS	ACNC	PT	URNS	ORD	MICROSCOPY.LIGHT	Mucous Threads UrnS Ql
5802-4	NITRITE	ACNC	PT	UR	ORD	TEST STRIP	Nitrite Ur Ql Strip
12453-7	PHOSPHATE CRYSTALS.AMORPHOUS	ACNC	PT	URNS	ORD	MICROSCOPY.LIGHT	Amorph Phos Cry UrnS Ql
2887-8	PROTEIN	ACNC	PT	UR	ORD		Prot Ur Ql
2888-6	PROTEIN	MCNC	PT	UR	QN		Prot Ur-mCnc
5812-3	SULFONAMIDE CRYSTALS	ACNC	PT	URNS	ORD	MICROSCOPY.LIGHT	Sulfonamide cry UrnS Ql
8249-5	TRANSITIONAL CELLS	ACNC	PT	URNS	ORD	MICROSCOPY.LIGHT	Trans Cells UrnS Ql
5813-1	TRICHOMONAS VAGINALIS	ACNC	PT	URNS	ORD	MICROSCOPY.LIGHT	T vaginalis UrnS Ql
5814-9	TRIPLE PHOSPHATE CRYSTALS	ACNC	PT	URNS	ORD	MICROSCOPY.LIGHT	Tri-Phos Cry UrnS Ql
5815-6	TYROSINE CRYSTALS	ACNC	PT	URNS	ORD	MICROSCOPY.LIGHT	Tyrosine Cry UrnS Ql
5817-2	URATE CRYSTALS	ACNC	PT	URNS	ORD	MICROSCOPY.LIGHT	Urate Cry UrnS Ql
12454-5	URATE CRYSTALS.AMORPHOUS	ACNC	PT	URNS	ORD	MICROSCOPY.LIGHT	Amorph Urate Cry UrnS Ql
5819-8	WAXY CASTS	NARIC	PT	URNS	QN	MICROSCOPY.LIGHT.LPF	Waxy Casts UrnS Qn LPF
21033-6	YEAST.BUDDING	ACNC	PT	URNS	ORD		Yeast Budding UrnS Ql

LOINC[®] Copyright Notice and License

The LOINC[®](R) codes, LOINC[®](R) Users' Guide, and LOINC[®](R) database are copyright (c) 1995-2005, Regenstrief Institute, Inc. and the Logical Observation Identifiers Names and Codes (LOINC[®]) Committee. All rights reserved.

LOINC[®](R) is a registered trademark of Regenstrief Institute, Inc.

Permission is hereby granted in perpetuity, without payment of license fees or royalties, to use, copy, or distribute the LOINC[®](R) codes, LOINC[®](R) Users' Guide, and the LOINC[®](R) database for any commercial or non-commercial purpose (RELMA(R) computer program, RELMA(R) database, and RELMA(R) Users' Guide, which are covered by a different license, are not included in this license), subject to the following terms and conditions:

1) To prevent the dilution of the purpose of the LOINC® codes and LOINC® database of providing a definitive standard for identifying clinical information in electronic reports, users shall not use the LOINC® Users' Guide, the LOINC® codes, and/or the LOINC® database for the purpose of developing or promulgating a different standard for identifying laboratory test results, diagnostic study reports, or clinical measurements and observations or orders for such entities in electronic reports and messages.

2) Users shall not change the meaning of any of the LOINC® codes. Users shall not change the name of, or any contents of, any fields in the LOINC® database. Users may add new fields to the LOINC® database to attach additional information to existing LOINC® records.

3) A user may delete records from the LOINC® database to deal with the user's local requirements. A user also may add new records to the LOINC® database to deal with the users' local requirements, provided that if new records are added, any new entry in the LOINC®_NUM field (field #1) of such new records must contain a leading alphabetic "X" so that the new codes and records cannot be confused with existing LOINC® codes or new LOINC® codes as they are defined in later releases of the LOINC® database. Records deleted or added by users to deal with local requirements are not reflected in the official LOINC® database maintained by the Regenstrief Institute and the LOINC® Committee. Users will also make reasonable efforts to submit requests to LOINC® for new records to cover observations that are not found in the LOINC® data base in order to minimize the need for X-codes.

4) LOINC® codes and other information from the LOINC® database may be used in electronic messages for laboratory test results and clinical observations such as HL7 ORU messages, without the need to include this LOINC® Copyright Notice and License or a reference thereto in the message (and without the need to include all fields required by Section 6 hereof). When the LOINC® code (from the LOINC®_NUM field) is included in the message, users are encouraged, but not required, to include the corresponding LOINC® short name (from the SHORTNAME field) in the message if the message provides a place for a text name representation of the code.

5) Users may make and distribute an unlimited number of copies of the unmodified LOINC® database. Each copy of the unmodified LOINC® database must include this LOINC® Copyright Notice and License, and must include the version number of the database. This LOINC® Copyright Notice and License must appear on every printed copy of the LOINC® database. Where the LOINC® database is distributed on a fixed storage medium (such as diskette or CD-ROM), a printed copy of this LOINC® Copyright Notice and License must be included on or with the storage medium, and a text file containing this information also must be stored on the storage medium in a file called "LOINC®_license.txt". Where the LOINC® database is distributed via the Internet, this LOINC® Copyright Notice and License must be accessible on the same Internet page from which the LOINC® database is available for download.

6) Subject to Section 1 and the other restrictions hereof, users may incorporate portions of the LOINC® database into another laboratory test definition database or software program for distribution outside of the user's corporation or organization, provided that any such laboratory test definition database or software program includes the following fields reproduced in their entirety from the LOINC® database: LOINC®_NUM (field #1), COMPONENT (field #2), PROPERTY (field #3), TIME_ASPCT (field #4), SYSTEM (field #5), SCALE_TYP (field #6), METHOD_TYP (field #7), ANSWERLIST (field #18), STATUS (field #19), and SHORTNAME (field #59). Users are encouraged, but not required, to also include the RelatedNames2 (field #58) in any such database. Further description of these fields is provided in Appendix A of the LOINC® Users' Guide. Every copy of the LOINC® database incorporated into or distributed in conjunction with another database or software program must include the following notice:

This product includes all or a portion of the LOINC®(R) database, or is derived from the LOINC®(R) database, subject to a license from Regenstrief Institute, Inc. Your use of the LOINC® database and LOINC® codes also is subject to this license, a copy of which is available at <http://www.LOINC.org/license>. The current complete LOINC® database and Users' Guide are available for download at <http://www.regenstrief.org/LOINC>. The LOINC® database and LOINC® codes are copyright (c) 1995-2005, Regenstrief Institute, Inc. and the Logical Observation Identifiers Names and Codes (LOINC®) Committee. All rights reserved. THE LOINC® DATABASE IS PROVIDED AS IS." ANY EXPRESS OR IMPLIED WARRANTIES ARE DISCLAIMED, INCLUDING, BUT NOT LIMITED TO, THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. LOINC®(R) is a registered trademark of Regenstrief Institute, Inc."

If the laboratory test definition database or software program containing the LOINC® database is distributed with a printed license, this statement must appear in the printed license. Where the laboratory test definition database or software program containing the LOINC® database is distributed on a fixed storage medium, a text file containing this information also must be stored on the storage medium in a file called "LOINC®_short_license.txt". Where the laboratory test definition database or software program containing the LOINC® database is distributed via the Internet, this information must be accessible on the same Internet page from which the product is available for download.

- 7) Use and distribution of the LOINC[®] codes and LOINC[®] database in ways that are not specifically discussed herein shall always be accompanied by the notice provided in Section 6 hereof. The guidelines for providing the notice that are contained in the last paragraph of Section 6 also shall apply. If a user has a question about whether a particular use of the LOINC[®] codes and/or the LOINC[®] database is permissible, the user is invited to contact the Regenstrief Institute by e-mail at LOINC[®]@regenstrief.org.
- 8) An unlimited number of copies of the LOINC[®] Users' Guide may be made and distributed. This LOINC[®] Copyright Notice and License must appear verbatim on every electronic or printed copy of the LOINC[®] Users' Guide. The LOINC[®] Users' Guide may not be modified, nor may derivative works of the LOINC[®] Users' Guide be created, without the prior written permission of the Regenstrief Institute, Inc. To request written permission, please contact LOINC[®]@regenstrief.org. The Regenstrief Institute retains the right to approve any modification to, or derivative work of, the LOINC[®] Users' Guide.
- 9) The names "Regenstrief," "Regenstrief Foundation," "Regenstrief Institute," and "LOINC[®] Committee" may not be used in a way which could be interpreted as an endorsement or a promotion of any product or service without prior written permission of the Regenstrief Institute, Inc. To request written permission, please contact LOINC[®]@regenstrief.org.
- 10) **DISCLAIMER: REGENSTRIEF INSTITUTE, INC. AND THE LOINC[®] COMMITTEE DO NOT ACCEPT LIABILITY FOR ANY OMISSIONS OR ERRORS IN THE LOINC[®] CODES, LOINC[®] USERS'GUIDE OR THE LOINC[®] DATABASE. THE LOINC[®] CODES, LOINC[®] USERS' GUIDE AND LOINC[®] DATABASE ARE PROVIDED "AS IS," WITHOUT WARRANTY OF ANY KIND. ANY EXPRESS OR IMPLIED WARRANTIES ARE DISCLAIMED, INCLUDING, BUT NOT LIMITED TO, THE IMPLIED WARRANTIES OF TITLE, NON-INFRINGEMENT, MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE OR WARRANTIES ARISING FROM A COURSE OF DEALING, TRADE USAGE,OR TRADE PRACTICE. FURTHER, NO WARRANTY OR REPRESENTATION IS MADE CONCERNING THE ACCURACY, COMPLETENESS, SEQUENCE, TIMELINESS OR AVAILABILITY OF THE LOINC[®] CODES, LOINC[®] USERS' GUIDE OR THE LOINC[®] DATABASE. IN NO EVENT SHALL REGENSTRIEF INSTITUTE, INC. OR THE LOINC[®] COMMITTEE OR ITS CONTRIBUTORS BE LIABLE FOR ANY DIRECT, INDIRECT, INCIDENTAL, SPECIAL, EXEMPLARY, RELIANCE, OR CONSEQUENTIAL DAMAGES OR ATTORNEYS' FEES (INCLUDING, BUT NOT LIMITED TO, PROCUREMENT OF SUBSTITUTE GOODS OR SERVICES; OPPORTUNITY COSTS; LOSS OF USE, DATA, SAVINGS OR PROFITS; OR BUSINESS INTERRUPTION) HOWEVER CAUSED AND ON ANY THEORY OF LIABILITY WHETHER IN CONTRACT, STRICT LIABILITY, OR TORT (INCLUDING NEGLIGENCE OR OTHERWISE) ARISING IN ANY WAY OUT OF THE USE OF LOINC[®] CODES, LOINC[®] USERS'GUIDE OR THE LOINC[®] DATABASE, EVEN IF ADVISED OF THE POSSIBILITY OF SUCH DAMAGE. SOME JURISDICTIONS DO NOT ALLOW THE LIMITATION OR EXCLUSION OF CERTAIN WARRANTIES OR CONDITIONS, SO SOME OF THE FOREGOING MAY NOT APPLY TO YOU.**
- 11) This license shall be construed and interpreted in accordance with the laws of the State of Indiana, United States of America, excluding its conflicts of law rules.