

# Indianapolis VA LOINC Mapping

Battery name	Test name	LOINC shortname	units	specimen	time	Data type	reflow	Help Text	LOINC
PLASMINOGEN ACTIVATOR INHIBITOR ACTIVITY	PLASMINOGEN ACTIVATOR INHIBITOR ACTIVITY			BLOOD	POINT	FREE TEXT		ANSWER MUST BE 1-15 CHARACTERS IN LENGTH	28165-9
MID %	MID %		%	BLOOD	POINT	FREE TEXT	0	ANSWER MUST BE 1-10 CHARACTERS IN LENGTH	28165-9
MISC DRUG SCREEN	OTHER DRUGS			URINE	POINT	SET		POS.P;NEG.N;TND.T; :B;	28165-9
AUTOIMMUNE LYMPHOCYTE TRANSFORMATION	AUTOIMMUNE LYMPHOCYTE TRANSFORMATION			BLOOD	POINT	FREE TEXT		ANSWER MUST BE 1-15 CHARACTERS IN LENGTH	28165-9
MISC DRUG SCREEN	OTHER DRUGS			BLOOD	POINT	SET	none det	POS.P;NEG.N;TND.T; :B;	28165-9
ZINC BLOOD	ZINC BLOOD	Zinc Bld-mCnc	ug/dL	BLOOD	POINT	NUMERIC	70	Number from 0 to 999 with 2 decimal	8245-3
XYLOSE	XYLOSE	Xylose Ur-mCnc		URINE	POINT	FREE TEXT		ANSWER MUST BE 1-7 CHARACTERS IN LENGTH	3135-1
XYLOSE	XYLOSE	Xylose Bld-mCnc		BLOOD	POINT	FREE TEXT		ANSWER MUST BE 1-7 CHARACTERS IN LENGTH	3134-4
SEMEN ANALYSIS	SEMEN WBC	WBC Srm Ql		SEMINAL FLUID	POINT	SET		FEW:F;MOD:M;MANY:Y;	32825-2
CBC	WBC	WBC # Ur Manual	#/uL	URINE	POINT	NUMERIC	0	Number from 0 to 999999 with 1 decimal	24122-4
CBC	WBC	WBC # Srv Manual	#/uL	SYNOVIAL FLUID	POINT	NUMERIC		Number from 0 to 999999 with 1 decimal	814-4
CBC	WBC	WBC # Pri Manual	#/uL	PERITONEAL FLUID	POINT	NUMERIC		Number from 0 to 999999 with 1 decimal	810-2
CBC	WBC	WBC # Pk Manual	#/uL	LUNG	POINT	NUMERIC		Number from 0 to 999999 with 1 decimal	808-6
CBC	WBC	WBC # Pk Manual	#/uL	PLEURAL FLUID	POINT	NUMERIC	0	Number from 0 to 999999 with 1 decimal	808-6
CBC	WBC	WBC # CSF Manual	#/cmm	CEREBROSPINAL FLUID	POINT	NUMERIC	0	Number from 0 to 999999 with 1 decimal	806-0
CBC	WBC	WBC # Bld Auto	K/uL	BLOOD	POINT	NUMERIC	4.5	Number from 0 to 999999 with 1 decimal	6680-2
COUMADIN	COUMADIN	Warfarin SerPl-mCnc	ug/mL	BLOOD	POINT	NUMERIC		K->X=X[(X>1000)](X<0)](X?.E 1".?N.N) X	4098-0
FATTY ACID PROFILE OF LIPIDS	TOTAL W6	w6 FA Ser-sCnc	umol/L	BLOOD	POINT	NUMERIC	3	Number from 0 to 9999 with 1 decimal	35178-1
FATTY ACID PROFILE OF LIPIDS	TOTAL W3	w3 FA Ser-sCnc	umol/L	BLOOD	POINT	NUMERIC	0.2	Number from 0 to 9999 with 1 decimal	35178-3
VMA PANEL	VMA RATIO	VMAcreat Ur-mRto	mcg/mg CR	URINE	POINT	NUMERIC	0	Number from 0 to 1000 with 1 decimal	3124-5
VMA	VMA	VMA 24h Ur-mRate	mg/24hr	URINE	POINT	NUMERIC	0	Number from 0 to 999 with 2 decimal	3122-9

Lab Section: CHEM

SATURATED C26:0	SATURATED C26:0	VLCFA C26:0 SerPl-mCnc	ug/ml	BLOOD	POINT	NUMERIC	0.15	K->X=(X>10)(X<0)(X?E!:" 4(N,N) X	33403-7
25 OH VITAMIN D3	25 OH VITAMIN D3	Vit D25 SerPl-mCnc	ng/mL	BLOOD	POINT	NUMERIC		Number from 0 to 999 with 2 decimal	19896-3
1,25-DIHYDROXYVIT D3	1,25-DIHYDROXYVIT D3	Vit D1,25 SerPl-mCnc	pg/mL	BLOOD	POINT	FREE TEXT		ANSWER MUST BE 1-7 CHARACTERS IN LENGTH	16496-3
VITAMIN C	VITAMIN C	Vit C Bld-mCnc	mg/dL	BLOOD	POINT	NUMERIC	0.4	S Q9="1,20,2" D "LRNUM	19026-6
VITAMIN B-6	VITAMIN B-6	Vit B6 SerPl-mCnc	ug/L	BLOOD	POINT	FREE TEXT		ANSWER MUST BE 1-6 CHARACTERS IN LENGTH	29006-9
VITAMIN B-2	VITAMIN B-2	Vit B2 Bld-mCnc	nmol/L	BLOOD	POINT	FREE TEXT	0.9	ANSWER MUST BE 1-6 CHARACTERS IN LENGTH	66896-1
VITAMIN B-12 (RIA)	VITAMIN B-12 (RIA)	Vit B12 Ser-mCnc	pg/mL	BLOOD	POINT	FREE TEXT	200	ANSWER MUST BE 1-6 CHARACTERS IN LENGTH	21326-9
VITAMIN B-12	VITAMIN B-12	Vit B12 Ser-mCnc	pg/mL	BLOOD	POINT	NUMERIC	239	Number from 0 to 100000 with 2 decimal	21326-9
VITAMIN B-1 (WHOLE BLOOD)	VITAMIN B-1 (WHOLE BLOOD)	Vit B1 Bld-sCnc	nmol/L	BLOOD	POINT	FREE TEXT		ANSWER MUST BE 1-6 CHARACTERS IN LENGTH	32554-8
VITAMIN B-1 (PLASMA)	VITAMIN B-1 (PLASMA)	Vit B1 Bld-mCnc	ug/L	BLOOD	POINT	FREE TEXT		ANSWER MUST BE 1-6 CHARACTERS IN LENGTH	29986-3
VITAMIN A	VITAMIN A	Vit A Bld-mCnc	ug/dL	BLOOD	POINT	NUMERIC		Number from 0 to 999 with 2 decimal	29226-3
SEMEN ANALYSIS	VISCOSITY	Viscosity Smm Ql		SEMINAL FLUID	POINT	FREE TEXT		ANSWER MUST BE 1-40 CHARACTERS IN LENGTH	32769-0
SEMEN ANALYSIS	VISCOSITY	Viscosity Ser Qn		SERUM	POINT	FREE TEXT	1.4	ANSWER MUST BE 1-40 CHARACTERS IN LENGTH	31286-6
SEMEN ANALYSIS	VISCOSITY	Viscosity Bld Qn		BLOOD	POINT	FREE TEXT	1.4	ANSWER MUST BE 1-40 CHARACTERS IN LENGTH	46804-4
VASOACTIVE PEPTIDE	VASOACTIVE PEPTIDE	VIP SerPl-mCnc	pg/mL	BLOOD	POINT	FREE TEXT	1	ANSWER MUST BE 1-6 CHARACTERS IN LENGTH	31256-2
CALCIUM CHANNEL BIND AB PQ-TYPE	CALCIUM CHANNEL BIND AB PQ-TYPE	VGCC-PQ Binding Ab Ser-sCnc		BLOOD	POINT	FREE TEXT		ANSWER MUST BE 1-10 CHARACTERS IN LENGTH	33980-4
CALCIUM CHANNEL BIND AB N-TYPE	CALCIUM CHANNEL BIND AB N- TYPE	VGCC-N Binding Ab Ser-sCnc		BLOOD	POINT	FREE TEXT		ANSWER MUST BE 1-10 CHARACTERS IN LENGTH	33979-6
VERAPAMIL	VERAPAMIL	Verapamil SerPl-mCnc	UG/L	BLOOD	POINT	FREE TEXT		ANSWER MUST BE 1-6 CHARACTERS IN LENGTH	40944-9
VANCOMYCIN,TROUGH	VANCOMYCIN,TROUGH	Vancomycin Tt SerPl-mCnc	mcg/mL	BLOOD	POINT	NUMERIC	5	Number from 0 to 9999 with 1 decimal	40926-3
FATTY ACID PROFILE OF LIPIDS	VACCENIC ACID	Vaccenate SerPl-sCnc	umol/L	BLOOD	POINT	NUMERIC	280	Number from 0 to 9999 with 1 decimal	35167-6
COMPUTED UREA NITROGEN	COMPUTED UREA NITROGEN	UUN 24H Ur-mRate	g/24HRS	URINE	POINT	FREE TEXT	6	Answer must be 1-7 characters in length.	30966-5
PORPHYRINS, FRACT.- FECES	UROPORPHYRIN	Uropor Ur Ql		URINE	POINT	FREE TEXT		ANSWER MUST BE 1-4 CHARACTERS IN LENGTH	31120-0
PORPHYRINS, FRACT.- FECES	UROPORPHYRIN	Uropor 24H Sit-mRate	ug/24 hr	FECES	POINT	FREE TEXT	< 80	ANSWER MUST BE 1-4 CHARACTERS IN LENGTH	17381-5

COMPUTED URIC ACID	COMPUTED URIC ACID	Urate 24H Ur-mRate	mg/24 hr	URINE	POINT	FREE TEXT	250	ANSWER MUST BE 1-7 CHARACTERS IN LENGTH	3087-4
URANIUM	URANIUM	Uranium Ur-mCnc		URINE	POINT	FREE TEXT		ANSWER MUST BE 1-10 CHARACTERS IN LENGTH	5752-1
TSH (STIMULATION)	TSH (STIMULATION)	TSH Recep Ab Ser-aCnc	miU/ml	BLOOD	POINT	FREE TEXT		ANSWER MUST BE 1-7 CHARACTERS IN LENGTH	5385-0
TRYPTASE	TRYPTASE	Tryptase SerPl-mCnc		BLOOD	POINT	FREE TEXT		ANSWER MUST BE 1-10 CHARACTERS IN LENGTH	21582-2
ELECSYS CARDIAC QC PANEL	TROPONIN-T	Troponin T SerPl-mCnc	ng/ml	PLASMA	POINT	NUMERIC	0	Number from 0 to 9999 with 2 decimal	6588-7
ELECSYS CARDIAC QC PANEL	TROPONIN-T	Troponin T BidV-mCnc	ng/ml	BLOOD	POINT	NUMERIC	0	Number from 0 to 9999 with 2 decimal	6587-9
TROPONIN-I	TROPONIN-I	Troponin I SerPl-mCnc	ng/ml	BLOOD	POINT	FREE TEXT	0	ANSWER MUST BE 1-10 CHARACTERS IN LENGTH	10839-9
TRI PHOSPHATE CRYSTALS	TRI PHOSPHATE CRYSTALS	Tri-Phos Cry Ums Ql		URINE	POINT	SET		NONE;NONE;RARE;RARE;OC C;OC;FEW;FEW;MOD;MOD; MANY;MANY;	5814-8
BETA-LOW DENSITY TRIGLYCERIDES	BETA-LOW DENSITY TRIGLYCERIDES	TrigI-VLDL Ester SerPl-mCnc		BLOOD	POINT	NUMERIC		Number from 0 to 9999 with 0 decimal	3047-8
LIPOPROTEIN PROFILE	Lp(e) Triglycerides	TrigI Bid-mCnc		BLOOD	POINT	FREE TEXT		ANSWER MUST BE 2-15 CHARACTERS IN LENGTH	3043-7
PHENOTHIAZINE SCREEN	TRIFLUOPROMAZINE	TriFluopromazine Ur Ql		URINE	POINT	FREE TEXT	neg	ANSWER MUST BE 1-7 CHARACTERS IN LENGTH	12322-4
PHENOTHIAZINE SCREEN	TRIFLUOPERAZINE	TriFluoperazine Ur Ql		URINE	POINT	FREE TEXT	neg	ANSWER MUST BE 1-7 CHARACTERS IN LENGTH	12321-6
COBAS BENZTCA QC PANEL	TRICYCLICS	Tricydics Ur Ql Scn		URINE	POINT	SET		NEG;NEG;POS;POS;	18312-8
ANTIDEPRESSANTS	ANTIDEPRESSANTS	Tricydics Ur Ql		URINE	POINT	SET		NEG;N;POS;P;SEE;S;TND;T; .B;	11004-9
ANTIDEPRESSANTS (SEND)	ANTIDEPRESSANTS (SEND)	Tricydics Ur Ql		URINE	POINT	SET		NEG;N;POS;P;TND;T;SEE;S; .B;	11004-9
COBAS BENZTCA QC PANEL	TRICYCLICS	Tricydics SerPl Ql		BLOOD	POINT	SET	none det	NEG;NEG;POS;POS;	4073-3
TRAZODONE	TRAZODONE	Trazodone SerPl-mCnc	ug/L	BLOOD	POINT	FREE TEXT	800	ANSWER MUST BE 1-6 CHARACTERS IN LENGTH	4064-2
NICOTINE PANEL	TRANS-3-OH-COTININE	trans-3-OH-cotinine Ur-mCnc	ng/ml	URINE	POINT	NUMERIC	<50	Number from 0 to 99999 with 1 decimal	33916-8
NICOTINE PANEL	TRANS-3-OH-COTININE	trans-3-OH-cotinine SerPl-mCnc	ng/ml	BLOOD	POINT	NUMERIC	<2.0	Number from 0 to 99999 with 1 decimal	35140-3
DIFF COUNT (BLOOD)	TOXIC GRANULATION	Toxic Granules Bid Ql Smear	sl/mkd	BLOOD	POINT	SET		Present;Present;	803-7
BRONCHIAL ALVEOLAR LAVAGE	CELL NUMBER (BAL)	Total Cells Counted Fld Qn	cells	BRONCHIAL WASHING CYTOLOGIC MATERIAL	POINT	FREE TEXT		ANSWER MUST BE 1-15 CHARACTERS IN LENGTH	38256-4

BRONCHIAL ALVEOLAR LAVAGE	CELLS/mL (BAL)	Total Cells Counted Fld Qn	cells/mL	BRONCHIAL WASHING CYTOLOGIC MATERIAL	POINT	FREE TEXT	ANSWER MUST BE 1-15 CHARACTERS IN LENGTH	38256-4
TOLUENE OCCUPATIONAL EXPOSURE	TOLUENE OCCUPATIONAL EXPOSURE	Toluene Bid-mCnc	mg/L	BLOOD	POINT	NUMERIC	Number from 0 to 9999 with 2 decimal	9637-0
TOCAINIDE	TOCAINIDE	Tocainide SerPl-mCnc	mcg/mL	BLOOD	POINT	NUMERIC	K*X=X(X>100)(X<0)(X?.E1" .3N) X	4068-0
TOBRAMYCIN,TROUGH	TOBRAMYCIN,TROUGH	Tobramycin Tr SerPl-mCnc	mcg/mL	BLOOD	POINT	NUMERIC	Number from 0 to 9999 with 1 decimal	4059-2
TRIDIONE	TRIDIONE	TMO SerPl-mCnc	mcg/dL	BLOOD	POINT	NUMERIC	K*X=X(X>100)(X<0)(X?.E1" .3N) X	4078-2
TIMENTIN	TIMENTIN	Ticarcillin SerPl-mCnc	UG/ML	BLOOD	POINT	FREE TEXT	Answer must be 1-8 characters in length.	4054-3
CHELATED IRON	CHELATED IRON	TIBC SerPl-mCnc	mcg/dL	BLOOD	POINT	FREE TEXT	ANSWER MUST BE 1-7 CHARACTERS IN LENGTH	2500-7
BETA-THROMBOGLOBULIN	BETA-THROMBOGLOBULIN	Thromboglob Pias-mCnc	mcg/dL	BLOOD	POINT	FREE TEXT	ANSWER MUST BE 1-10 CHARACTERS IN LENGTH	4686-2
THIOHIXENE	THIOHIXENE	Thiohixene SerPl-mCnc	NG/ML	BLOOD	POINT	FREE TEXT	Answer must be 1-10 characters in length.	6686-9
THIORIDAZINE	THIORIDAZINE	Thioridazine SerPl QI	mg/dL	BLOOD	POINT	FREE TEXT	ANSWER MUST BE 1-6 CHARACTERS IN LENGTH	13569-0
THIOCYANATE	THIOCYANATE	Thiocyanate Ur-mCnc	mg/dL	URINE	POINT	FREE TEXT	ANSWER MUST BE 1-6 CHARACTERS IN LENGTH	3004-9
THIOCYANATE	THIOCYANATE	Thiocyanate Bid-mCnc	mg/L	BLOOD	POINT	FREE TEXT	ANSWER MUST BE 1-6 CHARACTERS IN LENGTH	3001-5
THALLIUM	THALLIUM	Thallium Ur-mCnc	mg/L	URINE	POINT	FREE TEXT	ANSWER MUST BE 1-6 CHARACTERS IN LENGTH	5745-5
THALLIUM	THALLIUM	Thallium Bid-mCnc	mg/L	BLOOD	POINT	FREE TEXT	ANSWER MUST BE 1-6 CHARACTERS IN LENGTH	5743-0
ANABOLIC STEROIDS	T/E RATIO	Testost/Epilestost Ur-mRto	ng/dL	URINE	POINT	FREE TEXT	ANSWER MUST BE 1-15 CHARACTERS IN LENGTH	31030-0
BOUND TESTOSTERONE	BOUND TESTOSTERONE	Testost Bnd SerPl-mCnc	ng/dL	BLOOD	POINT	FREE TEXT	ANSWER MUST BE 1-7 CHARACTERS IN LENGTH	30123-4
FREE & WEAKLY BND TESTOST. %	FREE & WEAKLY BND TESTOST. %	Testost Bioavail % SerPl	%	BLOOD	POINT	FREE TEXT	ANSWER MUST BE 1-7 CHARACTERS IN LENGTH	6681-6
BENZODIAZEPINE SCREEN	TEMAZEPAM	Temazepam Ur QI Scn	neg	URINE	POINT	FREE TEXT	ANSWER MUST BE 1-7 CHARACTERS IN LENGTH	19688-0
BENZODIAZEPINE SCREEN	TEMAZEPAM	Temazepam SerPl QI	neg	BLOOD	POINT	FREE TEXT	ANSWER MUST BE 1-7 CHARACTERS IN LENGTH	13560-6
ACUTE LEUKEMIA PANEL	TERM.DEOXYNUCLEOTIDYL TRANS.	TdT Stn Bid/Mar		BLOOD	POINT	FREE TEXT	ANSWER MUST BE 1-20 CHARACTERS IN LENGTH	11021-3
ACUTE LEUKEMIA PANEL	TERM.DEOXYNUCLEOTIDYL TRANS.	TdT Stn Bid/Mar		BONE MARROW	POINT	FREE TEXT	ANSWER MUST BE 1-20 CHARACTERS IN LENGTH	11021-3
DIFF COUNT (BLOOD)	TARGET CELLS	Targets Bid QI Smear	s/mkd	BLOOD	POINT	SET	SLIGHT.SLT.MODERATE.MO D.MARKED.MKO.	10381-2

FK-506 LEVEL	FK-506 LEVEL	Tacrolimus Bid-mCnc	BLOOD	POINT	FREE TEXT	ANSWER MUST BE 1-10 CHARACTERS IN LENGTH	11253-2
THYROXIN BINDING GLOBULIN	THYROXIN BINDING GLOBULIN	T4BG Bid-mCnc	BLOOD	POINT	FREE TEXT	ANSWER MUST BE 1-6 CHARACTERS IN LENGTH	3020-5
FREE T4 (RIA)	FREE T4 (RIA)	T4 Free SerPI-mCnc	BLOOD	POINT	NUMERIC	Number from 0 to 9999 with 2 decimal	3024-7
FREE T-4 PANEL (SKBL-EQ.DIAL.)	FREE T-4 (EQ. DIAL.)	T4 Free SerPI Dialy-mCnc	BLOOD	POINT	FREE TEXT	Answer must be 1-8 characters in length.	6892-4
T UPTAKE	T UPTAKE	T3RU % SerPI	BLOOD	POINT	FREE TEXT	ANSWER MUST BE 1-5 NUMBERS	3050-2
REVERSE T-3	REVERSE T-3	T3Reverse SerPI-mCnc	BLOOD	POINT	NUMERIC	Number from 0 to 999 with 1 decimal	3052-8
FREE T3 (RIA)	FREE T3 (RIA)	T3Free SerPI-mCnc	BLOOD	POINT	NUMERIC	Number from 0 to 9999 with 2 decimal	3051-0
FREE T-3 INDEX	FREE T-3 INDEX	T3Free SerPI-mCnc	BLOOD	POINT	NUMERIC	Number from 1 to 9999 with 0 decimal	3051-0
FREE TRIIODOTHYRONINE (NEW)	FREE TRIIODOTHYRONINE (NEW)	T3Free SerPI-mCnc	BLOOD	POINT	FREE TEXT	ANSWER MUST BE 1-10 CHARACTERS IN LENGTH	3051-0
FREE T3 (MAYO)	FREE T3 (MAYO)	T3Free SerPI-mCnc	BLOOD	POINT	FREE TEXT	ANSWER MUST BE 1-15 CHARACTERS IN LENGTH	3051-0
TRIIODOTHYRONINE	TRIIODOTHYRONINE	T3 SerPI-mCnc	BLOOD	POINT	FREE TEXT	ANSWER MUST BE 1-5 NUMBERS	3053-6
TOTAL T-3	TOTAL T-3	T3 SerPI-mCnc	BLOOD	POINT	NUMERIC	Number from 0 to 9000 with 1 decimal	3053-6
FATTY ACID PROFILE OF LIPIDS	TRIENE TETRAENE RATIO	T:AA SerPI-eRto	BLOOD	POINT	NUMERIC	Number from 0 to 9999 with 3 decimal	35411-8
NEUROGENIC VIRAL BATTERY	TOXOPLASMA (IGM)	T gondii Igm Ser-eCnc	BLOOD	POINT	FREE TEXT	ANSWER MUST BE 1-8 CHARACTERS IN LENGTH	8040-8
NEUROGENIC VIRAL BATTERY	TOXOPLASMA (IGM)	T gondii Igm CSF Ql	CEREBROSPINAL FLUID	POINT	FREE TEXT	ANSWER MUST BE 1-8 CHARACTERS IN LENGTH	31181-1
NEUROGENIC VIRAL BATTERY	TOXOPLASMA (IGG)	T gondii Igg Titr Ser	BLOOD	POINT	FREE TEXT	ANSWER MUST BE 1-8 CHARACTERS IN LENGTH	22582-1
NEUROGENIC VIRAL BATTERY	TOXOPLASMA (IGG)	T gondii Igg CSF Ql	CEREBROSPINAL FLUID	POINT	FREE TEXT	ANSWER MUST BE 1-8 CHARACTERS IN LENGTH	31685-1
QUAL. DRUG ANALYSIS (SEND)	SYMPATHOMIMETIC AMINES (SEND)	Sympathomimetics Ur Ql	URINE	POINT	SET	NEG,N;POS;P;TND;T;SEE;S;.B;	29656-6
SULFONYLUREA HYPOGLY.-SCREEN	SULFONYLUREA HYPOGLY.-SCREEN	Sulfonylurea SerPI Ql	BLOOD	POINT	FREE TEXT	ANSWER MUST BE 1-10 CHARACTERS IN LENGTH	15374-2
SULFONAMIDES	SULFONAMIDES	Sulfonamide SerPI-mCnc	BLOOD	POINT	FREE TEXT	ANSWER MUST BE 1-6 CHARACTERS IN LENGTH	4040-2
SULFAPYRIDINE	SULFAPYRIDINE	Sulfapyridine SerPI-mCnc	BLOOD	POINT	FREE TEXT	ANSWER MUST BE 1-10 CHARACTERS IN LENGTH	14075-6
DIFF COUNT (BLOOD)	SPHEROCYTES	Spherocytes Bld Ql Smear	BLOOD	POINT	SET	RARE;RARE;FEW;FEW;SLIGH;T;S;T;MODERATE;MOD;MA;RKED;MKD;	802-9

SEMIN ANALYSIS	SPERM MORPHOLOGY	Spec Smn	SEMINAL FLUID	POINT	FREE TEXT	ANSWER MUST BE 2-58 CHARACTERS IN LENGTH	9704-8
SEMIN ANALYSIS	SPERM COUNT	Spec # Smn	SEMINAL FLUID	POINT	NUMERIC	Number from 0 to 999 with 0 decimal	9780-8
17-HYDROXYCORT. PANEL	VOLUME	Spec Vol Sfl	FECES	POINT	FREE TEXT	ANSWER MUST BE 1-8 CHARACTERS IN LENGTH	17606-5
17-HYDROXYCORT. PANEL	VOLUME	Spec Vol Smn	SEMINAL FLUID	POINT	FREE TEXT	ANSWER MUST BE 1-8 CHARACTERS IN LENGTH	3160-9
17-HYDROXYCORT. PANEL	VOLUME	Spec Vol Fid	GASTRIC CONTENTS	POINT	FREE TEXT	ANSWER MUST BE 1-8 CHARACTERS IN LENGTH	12254-9
17-HYDROXYCORT. PANEL	VOLUME	Spec Vol Dlaf	FLUID, DIALYSATE	POINT	FREE TEXT	ANSWER MUST BE 1-8 CHARACTERS IN LENGTH	12457-8
BRONCHIAL ALVEOLAR LAVAGE	VOLUME INSTILLED (BAL)	Spec Vol Bro	BRONCHIAL WASHING CYTOLOGIC MATERIAL	POINT	NUMERIC	Number from 1 to 9999 with 1 decimal	32801-3
BRONCHIAL ALVEOLAR LAVAGE	VOLUME RECIEVED (BAL)	Spec Vol Bro	BRONCHIAL WASHING CYTOLOGIC MATERIAL	POINT	NUMERIC	Number from 0 to 9999 with 1 decimal	32801-3
BRONCHIAL ALVEOLAR LAVAGE	VOLUME RETURNED (BAL)	Spec Vol Bro	BRONCHIAL WASHING CYTOLOGIC MATERIAL	POINT	NUMERIC	Number from 0 to 9999 with 1 decimal	32801-3
17-HYDROXYCORT. PANEL	VOLUME	Spec Vol 24H Ur-vRate	URINE	POINT	FREE TEXT	ANSWER MUST BE 1-8 CHARACTERS IN LENGTH	3168-2
BRONCHIAL ALVEOLAR LAVAGE	SOURCE (BAL)	Spec Source XXX	BRONCHIAL WASHING CYTOLOGIC MATERIAL	POINT	SET	LUL:Left Upper Lobe;LLL:Left Lower Lobe;RUL:Right Upper Lobe;RML:Right Middle Lobe;RLL:Right Lower Lobe;	31208-2
ANABOLIC STEROIDS	SPECIFIC GRAVITY	Sp Gr Ur Qn	URINE	POINT	NUMERIC	S X1=" S "<=>I\$E(X,1) X1=\$E(X,1);X=\$E(X,2,99) S:X=>X&(X,T,") X=X/1000+1;X=\$J(X,\$L(X),3) S X=X1_X S Q8="1,1.099,3"D ^LRNUM	2965-2
SPECIFIC GRAVITY (NON- URINE)	SPECIFIC GRAVITY (NON- URINE)	Sp Gr Srv Qn	SYNOVIAL FLUID	POINT	NUMERIC	S X1=" S "<=>I\$E(X,1) X1=\$E(X,1);X=\$E(X,2,99) S:X=>X&(X,T,") X=X/1000+1;X=\$J(X,\$L(X),3) S X=X1_X S Q8="1,1.099,3"D ^LRNUM	17598-4
SPECIFIC GRAVITY (NON- URINE)	SPECIFIC GRAVITY (NON- URINE)	Sp Gr Prt Qn Refractometry	PERTONEAL FLUID	POINT	NUMERIC	S X1=" S "<=>I\$E(X,1) X1=\$E(X,1);X=\$E(X,2,99) S:X=>X&(X,T,") X=X/1000+1;X=\$J(X,\$L(X),3) S X=X1_X S Q8="1,1.099,3"D ^LRNUM	14348-7

SPECIFIC GRAVITY (NON-URINE)	SPECIFIC GRAVITY (NON-URINE)	Sp Gr Pr Qn Refractometry	PLEURAL FLUID	POINT	NUMERIC	S X1=" S:"<=>TSE(X,1) X1=\$E(X,1);X=\$E(X,2,99) S X=>X&(X[".") X=X/1000+1;X=\$J(X,\$L(X),3) S X=X1_X S QB="1,1,099,3" D "LRNUM	14348-5
COMPUTED SODIUM	COMPUTED SODIUM	Sodium 24H Ur-sRate	URINE	POINT	FREE TEXT	ANSWER MUST BE 1-7 CHARACTERS IN LENGTH	2956-1
DIFF COUNT (BLOOD)	SICKLE CELLS	Sickle Cells Bld Ql Smear	BLOOD	POINT	SET	PRESENT;PRESENT;	901-1
SEX HORMONE BINDING GLOBULIN	SEX HORMONE BINDING GLOBULIN	SHBG Ser-mCnc	BLOOD	POINT	FREE TEXT	ANSWER MUST BE 1-10 CHARACTERS IN LENGTH	2942-1
SEROTONIN RELEASE	SEROTONIN RELEASE	Serotonin Bld-mCnc	BLOOD	POINT	FREE TEXT	ANSWER MUST BE 1-10 CHARACTERS IN LENGTH	2899-7
SEROTONIN-NEW	SEROTONIN-NEW	Serotonin Bld-mCnc	BLOOD	POINT	FREE TEXT	ANSWER MUST BE 1-10 CHARACTERS IN LENGTH	2899-7
SELENIUM	SELENIUM	Selenium Bld-mCnc	BLOOD	POINT	NUMERIC	Number from 0 to 999 with 2 decimal	5722-4
SECRETIN	SECRETIN	Secretin Plas-mCnc	BLOOD	POINT	FREE TEXT	ANSWER MUST BE 1-7 CHARACTERS IN LENGTH	2993-0
BARBITURATE SCREEN	SECOBARBITAL	Secobarbital Ur Ql	URINE	POINT	FREE TEXT	ANSWER MUST BE 1-8 CHARACTERS IN LENGTH	4028-5
BARBITURATE SCREEN	SECOBARBITAL	Secobarbital SerPl-mCnc	BLOOD	POINT	FREE TEXT	ANSWER MUST BE 1-8 CHARACTERS IN LENGTH	4028-7
DIFF COUNT (BLOOD)	SCHISTOCYTES	Schistocytes Bld Ql Smear	BLOOD	POINT	SET	RARE-RARE;FEW;FEW;SLUG HT;SLT;MODERATE;MOD;MA RKED;MKD;	800-3
FATTY ACID PROFILE OF LIPIDS	TOTAL SATURATED ACID	Sat FA SerPl-sCnc	BLOOD	POINT	NUMERIC	Number from 0 to 9999 with 1 decimal	35175-9
RAPID STREP	STREP SCREEN	S pyog Ag Thrt Ql	PHARYNX	POINT	FREE TEXT	ANSWER MUST BE 1-10 CHARACTERS IN LENGTH	18481-2
TOXIC SHOCK SYNDROME PANEL	TOXIC-SHOCK SYNDROME TOXIN-1	S aureus Toxt1 Ab Ser Ql	BLOOD	POINT	FREE TEXT	ANSWER MUST BE 1-10 CHARACTERS IN LENGTH	23936-8
NEUROGENIC VIRAL BATTERY	RESP SYNCYTIAL VIRUS ANTIGEN	RSV Ag XXX Ql	CEREBROSPINAL FLUID	POINT	FREE TEXT	ANSWER MUST BE 2-10 CHARACTERS IN LENGTH	31850-9
NEUROGENIC VIRAL BATTERY	RESP SYNCYTIAL VIRUS ANTIGEN	RSV Ab Titr Ser	BLOOD	POINT	FREE TEXT	ANSWER MUST BE 2-10 CHARACTERS IN LENGTH	22466-7
NEUROGENIC VIRAL BATTERY	TYPHUS AB TITER	Rick Typhus Ab Titr Ser	BLOOD	POINT	FREE TEXT	ANSWER MUST BE 2-10 CHARACTERS IN LENGTH	22482-3
SCAT	SCAT	RF Titr Ser	BLOOD	POINT	FREE TEXT	ANSWER MUST BE 1-7 CHARACTERS IN LENGTH	17857-4
RETINOL BINDING PROTEIN	RETINOL BINDING PROTEIN	Retinol Bind Prot SerPl-mCnc	BLOOD	POINT	FREE TEXT	ANSWER MUST BE 1-80 CHARACTERS IN LENGTH	1836-6
RENIN	RENIN	Renin Plas-cCnc	BLOOD	POINT	NUMERIC	Number from 0 to 999 with 2 decimal	2915-7
CBC	RDW	RDW RBC Auto-Rto	BLOOD	POINT	NUMERIC	Number from 0 to 100 with 1 decimal	788-0
SEMEN ANALYSIS	SEMEN RBC	RBC 5mn Ql Micro	SEMINAL FLUID	POINT	SET	FEW;F;MOD;M;MANY;Y;	13827-5

DIFF COUNT (BLOOD) CBC	MORPHOLOGY (RBC) RBC	RBC Bld	#/uL	BLOOD	POINT	SET	ESS NORM:ESS NORM; Number from 0 to 9999999 with 2 decimal	6742-1
CBC	RBC	RBC # Srv Manual	#/uL	SYNOVIAL FLUID	POINT	NUMERIC	Number from 0 to 9999999 with 2 decimal	797-1
CBC	RBC	RBC # Prt Manual	#/uL	PERITONEAL FLUID	POINT	NUMERIC	Number from 0 to 9999999 with 2 decimal	795-5
CBC	RBC	RBC # Plr Manual	#/uL	PLEURAL FLUID	POINT	NUMERIC	Number from 0 to 9999999 with 2 decimal	793-0
CBC	RBC	RBC # Fid Manual	#/uL	LUNG	POINT	NUMERIC	Number from 0 to 9999999 with 2 decimal	6741-3
CBC	RBC	RBC # CSF Manual	#/uL	CEREBROSPINAL FLUID	POINT	NUMERIC	Number from 0 to 9999999 with 2 decimal	792-2
CBC	RBC	RBC # Bld Auto	MuL	BLOOD	POINT	NUMERIC	Number from 0 to 9999999 with 2 decimal	789-8
NEUROGENIC VIRAL BATTERY	ROCKY MTN SPOTTED FV. TITER	Rickettsi Ab Titr Ser		BLOOD	POINT	FREE TEXT	ANSWER MUST BE 1-6 CHARACTERS IN LENGTH	29737-4
QUINIDINE	QUINIDINE	Quinidine Bld-mCnc	ug/mL	BLOOD	POINT	FREE TEXT	ANSWER MUST BE 0-7 CHARACTERS IN LENGTH	4014-7
PARANEOPLASTIC AUTOABY LUNG PANEL	PURKINJE CELL CYTOPLASMIC AB	Purkinje Cells Ab Ser Ql		BLOOD	POINT	FREE TEXT	ANSWER MUST BE 1-10 CHARACTERS IN LENGTH	27214-6
PTH (N-TERMINAL)	PTH (N-TERMINAL)	PTH-N-Term SerPl-mCnc	pg/mL	BLOOD	POINT	FREE TEXT	ANSWER MUST BE 1-7 CHARACTERS IN LENGTH	2733-4
PTH (INTACT)	PTH (INTACT)	PTH-Intact SerPl-mCnc	pg/ml	BLOOD	POINT	NUMERIC	Number from 0 to 5000 with 2 decimal	2731-8
PTH (C-TERMINAL)	PTH (C-TERMINAL)	PTH-C-Term SerPl-mCnc	pg/mL	BLOOD	POINT	FREE TEXT	ANSWER MUST BE 0-7 CHARACTERS IN LENGTH	2730-0
PTH (RELATED PROTEIN)	PTH (RELATED PROTEIN)	PTH Related Prot SerPl-sCnc	pmol/L	BLOOD	POINT	FREE TEXT	ANSWER MUST BE 1-7 CHARACTERS IN LENGTH	15087-0
PROTRIPTYLINE	PROTRIPTYLINE	Protip SerPl Ql		BLOOD	POINT	FREE TEXT	ANSWER MUST BE 1-6 CHARACTERS IN LENGTH	4000-6
PORPHYRINS, FRACT.- FECES	PROTOPORPHYRINS	Protopor RBC Ql		BLOOD	POINT	FREE TEXT	ANSWER MUST BE 1-10 CHARACTERS IN LENGTH	2897-7
PORPHYRINS, FRACT.- RBC	RBC PROTOPORPHYRIN	Protopor Bld-mCnc	MCG/Ml	BLOOD	POINT	FREE TEXT	Answer must be 1-8 characters in length.	2896-9
PORPHYRINS, FRACT.- FECES	PROTOPORPHYRINS	Protopor 24H Slt-mRate	mcg/24hrs	FECES	POINT	FREE TEXT	ANSWER MUST BE 1-10 CHARACTERS IN LENGTH	2891-0
CSF STUDIES	CSF STUDIES	Prot CSF-SerPl Ser+CSF	cells/mm3	CEREBROSPINAL FLUID	POINT	FREE TEXT	ANSWER MUST BE 1-60 CHARACTERS IN LENGTH	2858-9
CSF STUDIES	CSF STUDIES	Prot CSF Ql		BLOOD	POINT	FREE TEXT	ANSWER MUST BE 1-60 CHARACTERS IN LENGTH	2879-5
PROSTAGLANDIN E2, PLASMA	PROSTAGLANDIN E2, PLASMA	Prostaglandin E2 Ser-mCnc	PG/Ml	BLOOD	POINT	NUMERIC	Number from 0 to 9999 with 0 decimal	2853-0
PROPRANOLOL	PROPRANOLOL	Propranolol SerPl-mCnc	MCG/L	BLOOD	POINT	NUMERIC	Number from 0 to 2000 with 0 decimal	3982-5
PROPOXYPHENE CONFIRMATION	PROPOXYPHENE CONFIRMATION	Propoxyph Ur Ql Cfm		URINE	POINT	FREE TEXT	ANSWER MUST BE 1-10 CHARACTERS IN LENGTH	16200-8

DIFF COUNT (BLOOD)	PROMYELO	Promyelocytas % Bid Manual	%	BLOOD	POINT	NUMERIC	0	Number from 0 to 100 with 1 decimal	783-1
PROLACTIN (IU-2ND WHO)	PROLACTIN (IU-2ND WHO)	Prolactin SerPI-mCnc	ng/mL	BLOOD	POINT	FREE TEXT		ANSWER MUST BE 0-7 CHARACTERS IN LENGTH	2842-3
PROLACTIN	PROLACTIN	Prolactin SerPI-mCnc	ng/mL	BLOOD	POINT	FREE TEXT		ANSWER MUST BE 1-10 CHARACTERS IN LENGTH	2842-3
PROINSULIN	PROINSULIN	Proinsulin SerPI-sCnc	pmol/L	BLOOD	POINT	FREE TEXT	3	ANSWER MUST BE 1-10 CHARACTERS IN LENGTH	27882-0
PROGESTERONE	PROGESTERONE	Progester SerPI-mCnc	ng/mL	BLOOD	POINT	FREE TEXT		ANSWER MUST BE 1-5 CHARACTERS IN LENGTH	2839-9
ELECSYS CARDIAC QC PANEL	PRO-BNP	proBNP SerPI-mCnc	pg/mL	BLOOD	POINT	NUMERIC		Number from 0 to 99999 with 0 decimal	33762-6
PRIMIDONE	PRIMIDONE	Primidone SerPI-mCnc	mg/L	BLOOD	POINT	NUMERIC		Number from 0 to 999 with 2 decimal	3978-4
17-KETOSTEROID FRACTIONATION PANEL-UR	PREGNANETRIOL	Pregnanetriol 24H Ur-mRate	ug/24hr	URINE	POINT	FREE TEXT	245	ANSWER MUST BE 1-6 CHARACTERS IN LENGTH	2836-5
PREGNANEDIOL	PREGNANEDIOL	Pregnanediol 24H Ur-mRate	mg/24 hr	URINE	POINT	FREE TEXT		ANSWER MUST BE 1-6 CHARACTERS IN LENGTH	2834-0
LIPOPROTEIN PROFILE	PRE-BETA VERY LOW DEN. TRIG	PreB-LP SerPI QI		BLOOD	POINT	NUMERIC		Number from 0 to 9999 with 0 decimal	17847-5
LIPOPROTEIN PROFILE	PRE-BETA VERY LOW DEN. CHOL	PreB-LP SerPI QI		BLOOD	POINT	NUMERIC		Number from 0 to 9999 with 0 decimal	17847-5
LIPOPROTEIN ELECTROPHORESIS	LIPOPROTEIN, PRE-BETA	PreB-LP PIR QI	%	PLEURAL FLUID	POINT	FREE TEXT		ANSWER MUST BE 1-7 CHARACTERS IN LENGTH	27355-7
LIPOPROTEIN ELECTROPHORESIS	LIPOPROTEIN, PRE-BETA	PreB-LP % SerPI	%	SERUM	POINT	FREE TEXT	12.5	ANSWER MUST BE 1-7 CHARACTERS IN LENGTH	15123-3
PREALBUMIN	PREALBUMIN	Prealb SerPI-mCnc	mg/dl	BLOOD	POINT	FREE TEXT	20	ANSWER MUST BE 1-80 CHARACTERS IN LENGTH	14338-8
PHENYLPROPANOLAMINE	PHENYLPROPANOLAMINE	PPA SerPI QI		BLOOD	POINT	FREE TEXT		ANSWER MUST BE 1-7 CHARACTERS IN LENGTH	3963-6
COMPUTED POTASSIUM	COMPUTED POTASSIUM	Potassium 24H Ur-sRate	mmol/24 hr	URINE	POINT	FREE TEXT	40	ANSWER MUST BE 1-7 CHARACTERS IN LENGTH	2829-0
PORPHYRINS, QUANTITATIVE	PORPHYRINS, QUANTITATIVE	Porphyrins Ur-mCnc		URINE	POINT	FREE TEXT		ANSWER MUST BE 1-10 CHARACTERS IN LENGTH	11225-0
PORPHYRIN SCREEN	PORPHYRIN SCREEN	Porphyrins Ur QI		URINE	POINT	FREE TEXT		ANSWER MUST BE 1-8 CHARACTERS IN LENGTH	2818-3
PORPHYRINS, QUAN. PLASMA (PANEL)	PORPHYRINS, TOTAL PLASMA	Porphyrins SerPI-mCnc		BLOOD	POINT	NUMERIC		Number from 0 to 9999 with 1 decimal	2815-9
PORPHYRINS, QUANTITATIVE	PORPHYRINS, QUANTITATIVE	Porphyrins RBC-mCnc		BLOOD	POINT	FREE TEXT		ANSWER MUST BE 1-10 CHARACTERS IN LENGTH	2814-2
FATTY ACID PROFILE OF LIPIDS	TOTAL POLYUNSATURATED ACID	Polyunsat FA SerPI-sCnc	umol/L	BLOOD	POINT	NUMERIC	3.2	Number from 0 to 9999 with 1 decimal	35177-5
DIFF COUNT (BLOOD)	POLYCHROMASIA	Polychromasia Bid QI Smear	sl/mkd	BLOOD	POINT	SET		FEW;FEW;SLIGHT;SLT;MODE RATE;MOD;MARKED;MKD;	10378-8
DIFF COUNT (BLOOD)	POIKILOCYTOSIS	Poikilocytosis Bid QI Smear	sl/mkd	BLOOD	POINT	SET		COMMENT;COMMENT;	779-9

PO2 AT PT. TEMP	PO2 AT PT. TEMP	PO2 temp adj BldA Qn	ARTERIAL BLOOD	POINT	NUMERIC	Number from 15 to 600 with 1 decimal	19255-9
CBC	MPV	PMV Bld Qn Auto	BLOOD	POINT	NUMERIC	Number from 0 to 100 with 1 decimal	32823-1
PLASMINOGEN ASSAY	PLASMINOGEN ASSAY	PLG Act/Nor PPP Chro	BLOOD	POINT	FREE TEXT	ANSWER MUST BE 1-6 CHARACTERS IN LENGTH	28660-9
DIFF COUNT (BLOOD)	PLT (MORPH)	Platelet Bld	BLOOD	POINT	SET	ESS NORMLESS NORM/LARGE/LARGE/CLUMPED;	11125-2
DIFF COUNT (BLOOD)	PLT (ESTM)	Platelet Bld	BLOOD	POINT	SET	ADEQUATE/ADEQUATE/REDUCED/REDUCED/INCREASE/INC/HI ADEQ/LO ADEQ/LO ADEQ;	11125-2
CBC	PLT	Platelet # Bld Auto	BLOOD	POINT	NUMERIC	Number from 0 to 5000 with 1 decimal	777-3
DIFF COUNT (BLOOD)	PLASMA CELLS	Plasma Cells % Bld	BLOOD	POINT	NUMERIC	Number from 0 to 999 with 2 decimal	13047-6
PYRUVATE KINASE, ERYTHROCYTE	PYRUVATE KINASE, ERYTHROCYTE	PK RBC-sCnc	BLOOD	POINT	FREE TEXT	ANSWER MUST BE 1-15 CHARACTERS IN LENGTH	2912-4
PYRUVATE	PYRUVATE	PK Bld Ql	BLOOD	POINT	FREE TEXT	ANSWER MUST BE 1-8 CHARACTERS IN LENGTH	11227-6
PHYTTANATE	PHYTTANATE	Phytanate SerPl-mCnc	BLOOD	POINT	FREE TEXT	ANSWER MUST BE 1-10 CHARACTERS IN LENGTH	2804-3
COMPUTED PHOSPHORUS	COMPUTED PHOSPHORUS	Phosphate 24H Ur-mRate	URINE	POINT	FREE TEXT	ANSWER MUST BE 1-7 CHARACTERS IN LENGTH	2778-7
BARBITURATE SCREEN	PHENYTOIN	Phenytoin Ur Ql	URINE	POINT	FREE TEXT	ANSWER MUST BE 0-7 CHARACTERS IN LENGTH	12310-9
BARBITURATE SCREEN	PHENYTOIN	Phenytoin SerPl-mCnc	BLOOD	POINT	FREE TEXT	ANSWER MUST BE 0-7 CHARACTERS IN LENGTH	3988-5
FREE DILANTIN	FREE DILANTIN	Phenytoin Free SerPl-mCnc	BLOOD	POINT	FREE TEXT	ANSWER MUST BE 1-7 CHARACTERS IN LENGTH	3988-3
PHENOTHIAZINE SCREEN	PHENOTHIAZINES	Phenothiaz Ur Ql Scn	BLOOD	POINT	SET	POS:P;NEG:N;SEE:S;TND:T;.B;	19670-9
PHENOTHIAZINE SCREEN	PHENOTHIAZINES	Phenothiaz Ur Ql Scn	URINE	POINT	SET	POS:P;NEG:N;SEE:S;TND:T;.B;	19670-9
PHENOTHIAZINES (SEND)	PHENOTHIAZINES (SEND)	Phenothiaz Ur Ql	URINE	POINT	SET	NEG:N;POS:P;TND:T;SEE:S;.B;	3953-7
STOOL ALKALINIZATION	STOOL ALKALINIZATION	Phenolphthalein SH Ql	FECES	POINT	FREE TEXT	ANSWER MUST BE 3-15 CHARACTERS IN LENGTH	26868-0
PHENELZINE	PHENELZINE	Phenelzine SerPl-mCnc	BLOOD	POINT	FREE TEXT	ANSWER MUST BE 1-7 CHARACTERS IN LENGTH	9778-2
PHENACETIN	PHENACETIN	Phenac SerPl-mCnc	BLOOD	POINT	FREE TEXT	ANSWER MUST BE 1-7 CHARACTERS IN LENGTH	3929-7
ANABOLIC STEROIDS	PH	pH Ur-sCnc	URINE	POINT	FREE TEXT	ANSWER MUST BE 1-12 CHARACTERS IN LENGTH	2756-5
PH AT PT. TEMP	PH AT PT. TEMP	pH temp adj BldA-sCnc	BLOOD	POINT	NUMERIC	Number from 5 to 9 with 3 decimal	33254-4

ANABOLIC STEROIDS	PH	pH Srm-sCnc	SEMINAL FLUID	POINT	FREE TEXT	7.5	ANSWER MUST BE 1-12 CHARACTERS IN LENGTH	2752-4
ANABOLIC STEROIDS	PH	pH Prt-sCnc	PERTONEAL FLUID	POINT	FREE TEXT		ANSWER MUST BE 1-12 CHARACTERS IN LENGTH	2751-6
ANABOLIC STEROIDS	PH	pH Pfr-sCnc	PLEURAL FLUID	POINT	FREE TEXT		ANSWER MUST BE 1-12 CHARACTERS IN LENGTH	2750-8
ANABOLIC STEROIDS	PH	pH Gaat-sCnc	GASTRIC CONTENTS	POINT	FREE TEXT	5	ANSWER MUST BE 1-12 CHARACTERS IN LENGTH	2749-0
ANABOLIC STEROIDS	PH	pH Gaat-sCnc	PERICARDIAL FLUID	POINT	FREE TEXT		ANSWER MUST BE 1-12 CHARACTERS IN LENGTH	2749-0
ANABOLIC STEROIDS	PH	pH Fld-sCnc	FLUID, DIALYSATE	POINT	FREE TEXT	5	ANSWER MUST BE 1-12 CHARACTERS IN LENGTH	2748-2
ANABOLIC STEROIDS	PH	pH CSF-sCnc	CEREBROSPINAL FLUID	POINT	FREE TEXT		ANSWER MUST BE 1-12 CHARACTERS IN LENGTH	2747-4
ANABOLIC STEROIDS	PH	pH BidV-sCnc	BLOOD, VENOUS	POINT	FREE TEXT	7.35	ANSWER MUST BE 1-12 CHARACTERS IN LENGTH	2746-6
ANABOLIC STEROIDS	PH	pH Bid-sCnc	BLOOD	POINT	FREE TEXT	7.35	ANSWER MUST BE 1-12 CHARACTERS IN LENGTH	11558-4
ANABOLIC STEROIDS	PH	pH BidA-sCnc	ARTERIAL BLOOD	POINT	FREE TEXT	7.35	ANSWER MUST BE 1-12 CHARACTERS IN LENGTH	2744-1
PERPHENAZINE	PERPHENAZINE	Perphenazine SerPl Ql	BLOOD	POINT	FREE TEXT		ANSWER MUST BE 1-10 CHARACTERS IN LENGTH	26846-6
BARBITURATE SCREEN	PENTOBARBITAL (QUAL)	Pentobarb Ur Ql	URINE	POINT	FREE TEXT	neg	ANSWER MUST BE 1-6 CHARACTERS IN LENGTH	3925-5
PENTOBARBITAL (QUAN)	PENTOBARBITAL (QUAN)	Pentobarb SerPl Ql	BLOOD	POINT	NUMERIC		Number from 0 to 1000 with 1 decimal	3923-0
BARBITURATE SCREEN	PENTOBARBITAL (QUAL)	Pentobarb SerPl Ql	BLOOD	POINT	FREE TEXT		ANSWER MUST BE 1-6 CHARACTERS IN LENGTH	3923-0
PENTACARBOXYLPORPHYRIN RIN	PENTACARBOXYLPORPHYRIN	Pentis-CP Ur Ql	URINE	POINT	FREE TEXT		ANSWER MUST BE 1-10 CHARACTERS IN LENGTH	11222-7
FAT CELL TOXIC PROFILE	FAT CELL TOXIC PROFILE	Pentachlorophen Fat-mCnt	ADIPOSE TISSUE	POINT	FREE TEXT		ANSWER MUST BE 1-12 CHARACTERS IN LENGTH	20937-9
aPE-IgM-Indep	aPE-IgM-Indep	PE IgM Ser-aCnc	BLOOD	POINT	FREE TEXT	<=5	ANSWER MUST BE 1-10 CHARACTERS IN LENGTH	13077-3
aPE-IgM-Dep	aPE-IgM-Dep	PE IgM Ser-aCnc	BLOOD	POINT	FREE TEXT	<=5	ANSWER MUST BE 1-10 CHARACTERS IN LENGTH	13077-3
aPE-IgG-Indep	aPE-IgG-Indep	PE IgG Ser-aCnc	BLOOD	POINT	FREE TEXT	<=4	ANSWER MUST BE 1-10 CHARACTERS IN LENGTH	13076-5
aPE-IgG-Dep	aPE-IgG-Dep	PE IgG Ser-aCnc	BLOOD	POINT	FREE TEXT	<=4	ANSWER MUST BE 1-10 CHARACTERS IN LENGTH	13076-5
aPE-IgA-Dep	aPE-IgA-Dep	PE IgA Ser-aCnc	BLOOD	POINT	FREE TEXT	<=4	ANSWER MUST BE 1-10 CHARACTERS IN LENGTH	13078-1
aPE-IgA-Indep	aPE-IgA-Indep	PE IgA Ser-aCnc	BLOOD	POINT	FREE TEXT	<=4	ANSWER MUST BE 1-10 CHARACTERS IN LENGTH	13078-1
PHENCYCLIDINE CONFIRMATION	PHENCYCLIDINE CONFIRMATION	PCP Ur QlCfm	URINE	POINT	FREE TEXT		ANSWER MUST BE 1-10 CHARACTERS IN LENGTH	18392-1

